
ACTICLOT® Protein S

REF 843L

EC

REP

Obelis s.a
Bd. Général Wahis 53, 1030 Brussels, BELGIUM

INTENDED USE

ACTICLOT® Protein S is a functional clotting assay intended for the quantitative determination of Protein S activity in human plasma. The assay is for *in vitro* diagnostic use.

PRINCIPLE

The ACTICLOT Protein S assay is a clotting-based plasma assay. In the assay, dilutions of patient plasma are mixed with Protein S depleted plasma. A reagent that contains factor Xa, activated Protein C and phospholipids is added to the mixture. Following a 5-minute incubation period, calcium chloride is added to initiate clot formation. Under these conditions, the prolongation of the clotting time is directly proportional to the concentration of Protein S in the patient plasma.

BACKGROUND

DiScipio and colleagues first identified human Protein S (S for Seattle) in 1977.¹ Human Protein S is a single chain glycoprotein of 635 amino acids (69,000 kD) that circulates in plasma at a concentration of about 320 nM. The active Protein S gene has 15 exons and spans more than 80 kilobases on Chromosome 3.² Protein S is a vitamin K-dependant protein primarily synthesized in the liver but also in the endothelium and possibly in megakaryocytes. It is possible that all three sources of Protein S play a significant role in regulating haemostasis.

Protein S did not appear to be the zymogen precursor of a serine proteinase and no function was known until 1980 when Walker reported that the bovine Protein S functions as a cofactor for activated Protein C in the inactivation of activated factor V.³ The biologic importance of the regulation of coagulation by the pathway involving Protein C, Protein S and thrombomodulin has been substantiated by clinical studies. An association between thrombosis and partial deficiencies of Protein S has been documented in certain families.^{4,5,6} Protein S has been also found on the surface of endothelial cells and platelets.⁷ Approximately 60% of Protein S in human plasma is tightly bound to C4b-binding protein, therefore, plasma levels of Protein S may not completely reflect its bioavailability.⁸ This complex has no anticoagulant cofactor activity. It has been suggested that Protein S may have a role in the complement system by mediating the association of C4b-binding protein (C4bBP) with cell surfaces that have been damaged or activated to expose negatively charged phospholipids.⁹ This may represent another important regulatory mechanism for Protein S.

REAGENTS

- R1 4 vials ACTICLOT Activator Reagent (each vial contains 1 mL lyophilised material containing human activated Protein C, bovine factor Xa, and rabbit brain phospholipid).
- R2 4 vials Protein S Depleted Plasma (each vial contains 1 mL of lyophilised human plasma that has been depleted of Protein S by immunoadsorption on a column of immobilized goat antibody to human Protein S).
- R3 2 vials Sample Dilution Buffer, 10X Concentrate (2.5 mL concentrate, which upon dilution contain 0.2 M NaCl, 0.03 M HEPES and a sufficient quantity of Polybrene® † to neutralize up to 1.2 U/mL heparin in a plasma sample).
- R4 2 vials Protein S Control Plasma (each vial contains 0.5 mL of citrated lyophilised normal human plasma whose Protein S clotting activity level has been determined and is printed on the vial label).

This kit contains sufficient reagents for 40-160 test points, depending on whether performed by manual or automated method and the type of instrument utilized.

WARNING

Each donor unit used in the preparation of the activator reagent, depleted and control plasmas has been tested by an FDA cleared method and found to be non-reactive for Hepatitis B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus Type 1 (HIV-1). As no known test method can offer complete assurance that HBsAg, HCV, HIV-1 or other infectious agents are absent, this reagent should be handled as recommended for any potentially infectious human serum or blood specimen. Observe good laboratory hygiene practices.

MATERIALS REQUIRED BUT NOT PROVIDED

Automated photo-optical coagulation analyser, or mechanical clot timer, or stopwatch
0.025 M Calcium chloride solution
Type I deionised water
Incubator or water bath at 37°C
25 mL graduated cylinder
100-1000 µL variable volume pipette
Graph paper
Centrifuge (capable of 2,500 x g)

REAGENT RECONSTITUTION AND STABILITY

Unreconstituted reagents are stable until the expiration date indicated on the vial label, when stored at 2° - 8° C.

A. ACTICLOT Activator Reagent (R1)

Reconstitute with 1 mL Type I deionised water. Let stand at room temperature (18° - 25°C) for 20 minutes for complete dissolution. Store at 2° - 8°C and use within 4 hours. Reconstituted Activator Reagent may be frozen for:

	30 days
	-20°C

B. Protein S Depleted Plasma (R2)

Reconstitute with 1 mL Type I deionised water. Mix gently and let stand at room temperature for 20 minutes for complete dissolution. Store at 2° - 8°C and use within 4 hours. Reconstituted Depleted Plasma may be frozen for:

	30 days
	-20°C

C. Sample Dilution Buffer (R3)

Dilute contents of one vial to 25 mL with Type I deionised water. Store for:

	30 days
	2° - 8°C

D. Protein S Control Plasma (R4)

Reconstitute with 0.5 mL Type I deionised water. Mix gently and let stand at room temperature for 20 minutes for complete dissolution. Transfer to a plastic test tube. Reconstituted reagent is stable for:

	1 hour
	2° - 8°C

SPECIMEN COLLECTION AND PREPARATION

See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition", CLSI Document H21-A5, Vol. 28, No. 5, 2008¹⁰. Plasma collection should be performed as follows:

1. Collect 9 volumes of blood in 1 volume of 0.1 M trisodium citrate.
2. Centrifuge the blood sample at 2,500 x g for 15 minutes.
3. Plasma should be stored at 2° to 8°C and assayed within 2 hours. Alternatively, plasma may be stored for:

	1 month
	-20°C

4. Frozen plasma may be thawed once at 37°C, 30 minutes before use.

TRACEABILITY OF CONTROL MATERIALS

Information regarding traceability of control materials is available upon request from BioMedica Diagnostics.¹¹

QUALITY CONTROL

The Protein S Control plasma included in the kit has been evaluated against an internationally recognized reference preparation and may be tested as part of the laboratory quality control program for this assay. Control plasma should be included in the run whenever freshly reconstituted reagents are used. The Protein S level obtained for each control plasma should fall within the range specified by the vendor or laboratory. If the control sample fails to yield Protein S levels within the ranges specified, the run should be repeated. Contact BioMedica Diagnostics if repeated runs of controls fail to yield Protein S levels within acceptable limits.

ASSAY CALIBRATION

A new standard curve must be constructed each time the assay is performed. Pooled normal plasma (PNP) from at least 20 normal donors, which has been collected in the same manner as the plasma to be tested, should be used for the preparation of the assay standards. Commercially prepared reference plasma in which the Protein S level has been determined may also be used. Prepare the Protein S standards, control and patient plasma samples just before testing, as follows:

Standards and Samples – Recommended Dilutions

Protein S Standard Curve	PNP or Reference Plasma	Dilution Buffer
100%	100 µL	900 µL
50%	500 µL of 100% standard	500 µL
0%†	None	1000 µL
Patient and Control Plasma Samples	Patient and Control Plasma	Dilution Buffer
100%	50 µL	450 µL

† On some automated instruments; the lowest point will be a minimal dilution (approx. 12.5%) and not a "true" 0% standard.

Store diluted standards and samples on wet ice and test within 30 minutes.

ASSAY PROCEDURE

1. Transfer ACTICLOT Activator Reagent and Calcium chloride to the 37°C reagent wells of a clot timer instrument. Incubate for at least two minutes.
2. Add 0.1 mL Protein S Depleted Plasma and 0.1 mL test plasma dilution to a coagulation cuvette. Incubate for 2 minutes at 37°C.
3. Add 0.1 mL prewarmed ACTICLOT Activator Reagent to the cuvette. Incubate for exactly 5 minutes at 37°C.
4. Add 0.1 mL prewarmed Calcium Chloride. Start the timer and record the clotting time. Run each test sample and control in duplicate and average the duplicate results.


Note: To ensure accurate and reproducible results, use calibrated pipetting devices and observe the recommended procedures with particular emphasis on incubation times and temperatures.

CALIBRATION CURVE

Plot the percent Protein S activity of the calibration standards (x-axis) versus the mean clotting time (y-axis) on linear graph paper. Draw the line of best fit between the resulting points. For the patient and control samples, determine the percent Protein S activity by interpolating from the standard curve. In those patient samples with lupus anticoagulants or with abnormally high Protein S activity in which multiple dilutions were used, correct the Protein S level with the dilution factor. If the corrected Protein S levels still do not agree, the patient sample most likely has a very strong inhibitor for which an accurate activity level cannot be determined with this assay.

Sample Calibration Curve

The calibration curve shown below is for example only. Users must construct a new standard curve each time the assay is performed. Do not use the sample curve below.

LIMITATIONS AND INTERFERENCES

The presence of heparin greater than 1.2 U/mL in the patient sample, or the presence of lupus anticoagulants in the sample, may affect the assay results by prolonging the clotting time, giving an artificially high Protein S value. These samples should be tested at multiple dilutions to correct this effect.¹¹

Functional assays for Protein S can be influenced by the presence of factor V Leiden mutation in the patient, as this condition inhibits the Protein S enhanced inactivation of factor Va by activated Protein C.

This test should not be used as the sole means of diagnosis. Standard coagulation testing regimen should be followed.

EXPECTED VALUES

The value for Protein S is usually expressed as a relative percent compared to the value of a standard or pooled normal plasma. The expected normal range for Protein S is 55-160%.¹² However, each laboratory should determine its own normal range.¹³ Males have 10%-15% more plasma Protein S than females, and the concentration increases somewhat with age in both genders.⁹ A decrease in Protein S is associated with an increased incidence of thromboembolism.¹⁴ A decrease in Protein S activity does not necessarily indicate a decrease in Protein S plasma concentration. As Protein S (total) is present in plasma as free protein and as protein bound to C4bBP, only the free Protein S acts as a cofactor for activated Protein C. When the Protein S activity is decreased, it is important to establish the plasma levels of both free Protein S and C4bBP bound Protein S. Three types of congenital Protein S deficiency have been described, as summarized in Table 1.¹⁵

Table 1

DEFICIENCY TYPE	TOTAL PROTEIN S	FREE PROTEIN S	PROTEIN S ACTIVITY
I	∇ ^b	∇	∇
II	↔ ^a	↔	∇
III	↔	∇	∇

a: ↔ = is normal

b: ∇ = is at a reduced level

Levels of Protein S may also be decreased in liver disease and during anticoagulant therapy. In addition, Protein S activity and free Protein S antigen may be reduced in inflammatory disease, in which the levels of C4bBP are elevated.

PERFORMANCE CHARACTERISTICS

The end-user should establish product performance characteristics for the specific method or instrumentation used. Studies using an automated analyser with photo-optical detection produced the results shown below¹¹.

Accuracy

A study performed using another commercially available Protein S kit gave a correlation coefficient of 0.868 for 78 samples tested, with a regression equation of $y = 0.913x + 17.2$.

Precision

The following coefficients of variation (CV's) were observed:

Protein S Level	Intra-Assay CV	Inter-Assay CV
Normal Plasma (87%)	3.0%	5.9%
Abnormal Plasma (42%)	2.3%	5.1%

Sensitivity

ACTICLOT Protein S is designed to give a linear standard curve for Protein S levels between 0 and 100% when plotted as directed.

[†]Polybrene® is a registered trademark of Abbott Laboratories.

REFERENCES

1. Discipio, R. G., *et al.* A comparison of human Prothrombin, factor IX (Christmas Factor) factor X (Stuart Factor) and Protein S. *Biochemistry* 1977, **16**: 698-706.
2. Schmidel, D. K. Organization of the human Protein S genes. *Biochemistry* 1990, **16**: 7845-7852.
3. Walker, F. J. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. *Journal of Biological Chemistry* 1980, **255**: 5521-5524.
4. Comp, P. C., *et al.* Familial Protein S deficiency is associated with recurrent thrombosis. *Journal of Clinical Investigations* 1984, **74**: 2082-2088.
5. Broekmans, A. W., *et al.* Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. *Thrombosis and Haemostasis* 1985, **53**: 273-277.

6. Kamiya, T., *et al.* Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
7. Schwarz, H. P., *et al.* Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**: 1452-1455.
8. Dahlbäck, Björn. Interaction between vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in Thrombosis and Hemostasis* 1984, **10**: 139-148.
9. Broze, G. J., *et al.* Biochemistry and Physiology of Protein C, Protein S, and Thrombomodulin. In: Colman, *et al.*, eds. *Hemostasis and Thrombosis: Basic Principles and Clinical Practice, Third Edition* 1994. Philadelphia, PA: J. B. Lippincott Company, 259-275.
10. *Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition*, CLSI Document H21-A5, Vol. 28, No. 5, 2008.
11. Data on file at Sekiusi Diagnostics, LLC, Stamford, CT 06902.
12. Boyer-Neumann, C., *et al.* Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
13. How to Define and Determine Reference Intervals in the Clinical Laboratory – second edition. *NCCLS/CLSI Document C28-A2* June 2000; Vol. 20, No. 13
14. Bertina, R. M. Hereditary Protein S deficiency. *Haemostasis* 1995, **15**: 241-246.
15. Hathaway, W. E. and Goodnight, Jr., S. H. *Disorders of Hemostasis and Thrombosis: A Clinical Guide*. McGraw-Hill 1993, p. 350.

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LITERATURHINWEISE

- Discipio, R. G., et al. A comparison of human Prothrombin, factor IX (Christmas Factor) factor X (Stuart Factor) and Protein S. *Biochemistry* 1977, **16**: 698-706.
- Schmidel, D. K. Organization of the human Protein S genes. *Biochemistry* 1990, **16**: 7845-7852.
- Walker, F.J., 1980. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. *Journal of Biological Chemistry*, 255, 5521-5524.
- Comp, P. C., et al. Familial Protein S deficiency is associated with recurrent thrombosis. *Journal of Clinical Investigations* 1984, **74**: 2082-2088.
- Broekmans, A. W., et al. Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. *Thrombosis and Haemostasis* 1985, **53**: 273-277.
- Kamiya, T., et al. Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
- Schwarz, H. P., et al. Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**: 1452-1455.
- Dahlbäck, Björn. Interaction between vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in Thrombosis and Hemostasis* 1984, **10**: 139-148.
- Broze, G. J., et al. Biochemistry and Physiology of Protein C, Protein S, and Thrombomodulin. In: Colman, et al., eds. *Hemostasis and Thrombosis: Basic Principles and Clinical Practice, Third Edition* 1994. Philadelphia, PA: J. B. Lippincott Company, 259-275.
- Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; *Approved Guideline-Fifth Edition*, CLSI Document H21-A5, Vol. 28, No. 5, 2008.
- Daten erhältlich von BioMedica Diagnostics, Stamford, CT 06902.
- Boyer-Neumann, C., et al. Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
- How to Define and Determine Reference Intervals in the Clinical Laboratory – second edition. *NCCLS/CLSI Document C28-A2* June 2000; Vol. 20, No. 13.
- Bertina, R. M. Hereditary Protein S deficiency. *Haemostasis* 1995, **15**: 241-246.
- Hathaway, W. E. and Goodnight, Jr., S. H. *Disorders of Hemostasis and Thrombosis: A Clinical Guide*. McGraw-Hill 1993, p. 350.

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REFERENCES

- 1. Discipio, R.G., et al. A comparison of human Prothrombin, factor IX (Christmas factor), factor X (Stuart Factor) and Protein S. Biochemistry 1977, **16**: 696-706.
- 2. Schmidel, D.K. Organization of the human Protein S genes. Biochemistry 1990, **16**: 7845-7852.
- 3. Walker, F.J. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. Journal of Biological Chemistry 1980, **255**: 5521-5524.
- 4. Comp, P.C., et al. Familial Protein S deficiency is associated with a recurrent thrombosis. Journal of Clinical investigations 1984, **74**: 2082-2088.
- 5. Broekmans, A. W., et al. Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. Thrombosis and Haemostasis 1985, **53**: 273-277.

6. Kamiya, T., et al Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
7. Schwarz, H. P., et al. Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**:1452-1455.
8. Dahlbäck, Björn Interaction between Vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in thrombosis and hemostasis* 1984, **10**:139-148.
9. Broze, G.J., et al. Biochemistry and Physiology of Protein C, Protein S and thrombomodulin. In: *Colman, et al., eds Hemostasis and Thrombosis: Basic principles and Clinical Practice, third Edition* 1994 Philadelphia, PA, J.B. Lippincott Company, 259-275.
10. *Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition*, CLSI Document H21-A5, Vol 28, No 5, 2008.
11. Data on file at Sekiusi Diagnostics, LLC, Stamford, CT 06902.
12. Boyer-Neumann, C., et al. Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
13. How to define and Determine reference intervals in the Clinical Laboratory - second edition NCCLS/CLSI Document C28-A2 June 2000, Vol 20, No 13.
14. Bertina, R.M. Hereditary Protein S deficiency. *Haemostasis* 1985; **15**: 241-245.
15. Hathaway, W.E. and Goodnight, Jr, S.H. *Disorders of haemostasis and thrombosis: A clinical Guide*. McGraw- Hill 1993, p 350.

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1. Discipio, R. G., *et al.* A comparison of human Prothrombin, factor IX (Christmas Factor) factor X (Stuart Factor) and Protein S. *Biochemistry* 1977, **16**: 698-706.
2. Schmidel, D. K. Organization of the human Protein S genes. *Biochemistry* 1990, **16**: 7845-7852.
3. Walker, F. J. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. *Journal of Biological Chemistry* 1980, **255**: 5521-5524.
4. Comp, P. C., *et al.* Familial Protein S deficiency is associated with recurrent thrombosis. *Journal of Clinical Investigations* 1984, **74**: 2082-2088.
5. Broekmans, A. W., *et al.* Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. *Thrombosis and Haemostasis* 1985, **53**: 273-277.
6. Kamiya, T., *et al.* Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
7. Sschwarz, H. P., *et al.* Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**: 1452-1455.
8. Dahlbäck, Björn. Interaction between vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in Thrombosis and Hemostasis* 1984, **10**: 139-148.
9. Broze, G. J., *et al.* Biochemistry and Physiology of Protein C, Protein S, and Thrombomodulin. In: Colman, *et al.*, eds. *Hemostasis and Thrombosis: Basic Principles and Clinical Practice, Third Edition* 1994. Philadelphia, PA: J. B. Lippincott Company, 259-275.
10. *Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines-Fifth Edition*, CLSI Document H21-A5, Vol. 28, No. 5, 2008.
11. Datos registrados en BioMedica Diagnostics, Stamford, CT 06902.
12. Boyer-Neumann, C., *et al.* Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
13. How to Define and Determine Reference Intervals in the Clinical Laboratory – second edition. *NCCLS/CLSI Document C28-A2* June 2000; Vol. 20, No. 13
14. Bertina, R. M. Hereditary Protein S deficiency. *Haemostasis* 1995, **15**: 241-246.
15. Hathaway, W. E. and Goodnight, Jr., S. H. *Disorders of Hemostasis and Thrombosis: A Clinical Guide*. McGraw-Hill 1993, p. 350.

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BIBLIOGRAFIA

- Discipio, R. G., et al. A comparison of human Prothrombin, factor IX (Christmas Factor) factor X (Stuart Factor) and Protein S. *Biochemistry* 1977, **16**: 698-706.
- Schmidel, D. K. Organization of the human Protein S genes. *Biochemistry* 1990, **16**: 7845-7852.
- Walker, F. J. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. *Journal of Biological Chemistry* 1980, **255**: 5521-5524.
- Comp, P. C., et al. Familial Protein S deficiency is associated with recurrent thrombosis. *Journal of Clinical Investigations* 1984, **74**: 2082-2088.
- Broekmans, A. W., et al. Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. *Thrombosis and Haemostasis* 1985, **53**: 273-277.
- Kamiya, T., et al. Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
- Schwarz, H. P., et al. Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**: 1452-1455.
- Dahlbäck, Björn. Interaction between vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in Thrombosis and Hemostasis* 1984, **10**: 139-148.
- Broze, G. J., et al. Biochemistry and Physiology of Protein C, Protein S, and Thrombomodulin. In: Colman, et al., eds. *Hemostasis and Thrombosis: Basic Principles and Clinical Practice, Third Edition* 1994. Philadelphia, PA: J. B. Lippincott Company, 259-275.
- Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition, CLSI Document H21-A5, Vol. 28, No. 5, 2008.
- Dati su file disponibili presso la BioMedica Diagnostics.
- Boyer-Neumann, C., et al. Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
- How to Define and Determine Reference Intervals in the Clinical Laboratory – second edition. *NCCLS/CLSI Document C28-A2* June 2000; Vol. 20, No. 13.
- Bertina, R. M. Hereditary Protein S deficiency. *Haemostasis* 1995, **15**: 241-246.
- Hathaway, W. E. and Goodnight, Jr., S. H. *Disorders of Hemostasis and Thrombosis: A Clinical Guide*. McGraw-Hill 1993, p. 350.

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BIBLIOGRAFIA

1. Discipio, R. G., *et al.* A comparison of human Prothrombin, factor IX (Christmas Factor) factor X (Stuart Factor) and Protein S. *Biochemistry* 1977, **16**: 698-706.
2. Schmidel, D. K. Organization of the human Protein S genes. *Biochemistry* 1990, **16**: 7845-7852.
3. Walker, F. J. Regulation of activated Protein C by a new protein. A possible function for bovine Protein S. *Journal of Biological Chemistry* 1980, **255**: 5521-5524.
4. Comp, P. C., *et al.* Familial Protein S deficiency is associated with recurrent thrombosis. *Journal of Clinical Investigations* 1984, **74**: 2082-2088.
5. Broekmans, A. W., *et al.* Hereditary Protein S deficiency and venous thromboembolism. A study in three Dutch families. *Thrombosis and Haemostasis* 1985, **53**: 273-277.
6. Kamiya, T., *et al.* Inherited deficiency of Protein S in a Japanese family with recurrent venous thromboembolism. A study of three generations. *Blood* 1986, **67**: 406-410.
7. Schwarz, H. P., *et al.* Identification and quantification of Protein S in human platelets. *Blood* 1985, **66**: 1452-1455.
8. Dahlbäck, Björn. Interaction between vitamin K-dependant Protein S and the complement C4b-binding protein. *Seminars in Thrombosis and Hemostasis* 1984, **10**: 139-148.
9. Broze, G. J., *et al.* Biochemistry and Physiology of Protein C, Protein S, and Thrombomodulin. In: Colman, *et al.*, eds. *Hemostasis and Thrombosis: Basic Principles and Clinical Practice, Third Edition* 1994. Philadelphia, PA: J. B. Lippincott Company, 259-275.
10. *Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition*, CLSI Document H21-A5, Vol. 28, No. 5, 2008.
11. Datos registrados en BioMedica Diagnostics, Stamford, CT 06902.
12. Boyer-Neumann, C., *et al.* Comparison of functional assays for Protein S. *Thrombosis and Haemostasis* 1993, **70**: 946-950.
13. How to Define and Determine Reference Intervals in the Clinical Laboratory – Second edition. *NCCLS/CLSI Document C28-A2* June 2000; Vol. 20, No. 13
14. Bertina, R. M. Hereditary Protein S deficiency. *Haemostasis* 1995, **15**: 241-246.
15. Hathaway, W. E. and Goodnight, Jr., S. H. *Disorders of Hemostasis and Thrombosis: A Clinical Guide*. McGraw-Hill 1993, p. 350.

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