
IMUBIND® Plasma PAI-1 ELISA

REF 822

Obelis s.a
 Bd. Général Wahis 53, 1030 Brussels, BELGIUM

INTENDED USE

The IMUBIND® Plasma PAI-1 ELISA is an enzyme-linked immunosorbent assay for the quantitative measurement of human Plasminogen Activator Inhibitor Type-1 (PAI-1) antigen in plasma. This assay is for *in vitro* diagnostic use.

EXPLANATION OF THE TEST

Plasminogen Activator Inhibitor Type-1 (PAI-1), a primary regulator of fibrinolysis, has been found in a number of different tissues, and cell types including macrophages/monocytes, hepatocytes, vascular endothelia, adipose tissue of the heart and lungs, and in platelets.^{1,2} The clinical interest in measuring PAI-1 in plasma is based on case studies in which levels of this serine protease inhibitor are associated with various thrombotic and fibrinolytic complications. Deficiency of PAI-1 activity is associated with bleeding disorders wherein the routine haemostatic screening tests are normal.^{3,4} High levels of PAI-1 activity are found in patients suffering from myocardial infarction, haemolytic uremic syndrome, and stroke.⁵⁻⁹ Levels of PAI-1 in the plasma of pregnant women are also correlated with gestational diabetes, reduced placental blood flow and preeclampsia.^{10,11,12} Patients with cirrhosis may also have elevated levels of PAI-1.¹³

PRINCIPLE OF THE PROCEDURE

The IMUBIND Plasma PAI-1 ELISA uses microwells coated with a mouse monoclonal antibody directed against human PAI-1. PAI-1 standards and plasma samples are added to the wells and the antibody captures the PAI-1 antigen present during an incubation period. A peroxidase conjugated antibody against human PAI-1 is added to the well and this conjugate binds to the captured PAI-1 molecules. All unbound material is washed away and a peroxidase reactive enzyme substrate, ortho-phenylenediamine (OPD) is added to the wells. The subsequent peroxidase/substrate reaction yields a yellow colored solution. Addition of sulfuric acid stops the reaction and turns the color of the solution color orange. The absorbance of the solution is measured at 490 nm. The absorbance is directly proportional to the amount of PAI-1 present in the sample.

REAGENTS

- R1 96 Monoclonal Antibody Anti-Human PAI-1 Coated Microwells with 1 Acetate Sheet to cover the microwells during assays
- R2 1 vial PET Buffer (PBS-EDTA-Tween 20, powder), 15.15 g
- R3 1 vial PAI-1 Standard, 50 ng/mL (lyophilized)
- R4 1 vial PAI-1 Standard, 0 ng/mL (lyophilized)
- R5 1 vial Detection Antibody, 200 µL (concentrate)
- R6 2 tablets OPD Substrate (each tablet contains 10 mg of OPD-HCl)

*There are sufficient reagents to assay 44 plasma samples and generate a four-point standard curve (both tested in duplicate), for a total of 96 tests. Samples may be patient, control or reference plasmas.

WARNINGS and PRECAUTIONS

The source material for reagents in this kit have been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using US Food and Drug Administration approved methods. As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

Not for internal use in humans or animals. Do not use the kit components beyond the stated expiration date. Do not mix reagents from different kits. Avoid microbial contamination of the reagents. Do not smoke, eat or drink in areas in which specimens or kit reagents are handled. Do not pipette reagents by mouth. Wear laboratory coat and disposable gloves throughout the test procedure and wash hands thoroughly afterwards. Avoid splashing or aerosol formation. Do not expose the OPD reagent to strong sunlight. Store OPD in the dark and do not allow it to come into contact with any oxidizing agents or metal parts.

PET buffer	Warning		H319, P264, P280, P305+P351+P338, P337+P313
Detection antibody	Warning		H319, P264, P280, P305+P351+P338, P337+P313
OPD substrate	Warning	 	CONT o-Phenylenediamine dihydrochloride $C_6H_4(NH_2)_2 \cdot 2HCl$ H302, H319, H317, H341, H351, H410, P264, P280, P202, P261, P273, P301+P312, P330, P305+P351+P338, P302+P352

Hazard Statements¹⁴:

- H302 Harmful if swallowed.
- H317 May cause an allergic skin reaction.
- H319 Causes serious eye irritation.
- H341 Suspected of causing genetic defects.
- H351 Suspected of causing cancer.
- H410 Very toxic to aquatic life with long lasting effects.

Precautionary Statements:	P202	Do not handle until all safety precautions have been read and understood.
	P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
	P264	Wash thoroughly after handling.
	P273	Avoid release to the environment.
	P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
	P301 + P312	IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
	P330	Rinse mouth.
	P302 + P352	IF ON SKIN: Wash with plenty of water.
	P305 + P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	P337 + P313	If eye irritation persists: Get medical advice/attention.

MATERIALS REQUIRED BUT NOT PROVIDED

0.22 micron filtered deionised water
 8-channel pipette covering 50 - 200 µL
 1-channel pipettes covering 10 - 1000 µL
 repetitive pipette 50 µL
 squeeze bottle, towels or thin sponge
 50 mL containers, e.g., plastic tubes
 microwell plate washer, microwell plate orbital shaker
 microwell plate reader set at a wavelength of 490 nm
 concentrated sulfuric acid (H₂SO₄)
 30% hydrogen peroxide (H₂O₂)
 centrifuge (capable of 1,500 x g)

REAGENT PREPARATION AND STORAGE

A. PET Buffer (R2)

Dissolve the contents of the PET Buffer vial in 1 L of filtered deionised water. Mix for 15 minutes. Prepared buffer may be stored at:

	1 month
	2° - 8°C

B. PAI-1 Standards (R3, R4)

- Add 0.5 mL filtered deionised water to the 0 ng/mL PAI-1 Standard and the 50 ng/mL PAI-1 Standard vials. Agitate gently for three minutes. Treat the standards as you would a plasma sample (see SPECIMEN COLLECTION AND PREPARATION). Standards should be used within one hour after reconstitution. Reconstituted standards may also be aliquotted into tightly capped tubes and stored at:

	6 months
	-70°C or colder

- Prepare a 12.5 ng/mL standard and a 25 ng/mL standard for the calibration curve by adding the 50 ng/mL PAI-1 Standard and the 0 ng/mL PAI-1 Standard in the following ratios. Mix each tube for a few seconds. Pipette the 50 ng/mL and 0 ng/mL standards directly from their vials.

PAI-1 Standard	Volume of 0 ng/mL PAI-1 Standard	Volume of 50 ng/mL PAI-1 Standard
12.5 ng/mL	75 µL	25 µL
25 ng/mL	50 µL	50 µL

C. Detection Antibody (R5)

Prepare working strength Detection Antibody by diluting 1:50 with PET Buffer. For running all 96 microwells at one time, dilute 120 µL of Detection Antibody up to 6 mL in PET Buffer. If all 96 microwells are not to be used, dilute 20 µL of Detection Antibody up to 1 mL in PET Buffer for each 16 microwells that will be used. Working strength Detection Antibody may be stored in the dark for:

	1 month
	-20°C

D. Substrate (R6)

Add 1 tablet OPD substrate to 3 mL of filtered deionised water. After the tablet has fully dissolved, add 21 mL of filtered deionised water to the dissolved tablet for a final volume of 24 mL. Each 4 mL of substrate solution is sufficient for 16 sample wells. Prepared substrate may be aliquotted and stored at:

	1 month
	-20°C

E. Stop Solution

Prepare a 4.5 M H₂SO₄ solution by adding 5 mL of concentrated sulfuric acid (95-97%) to 15 mL of filtered deionised water. Mix well. Store unused Stop Solution at:

	2° - 8°C
---	----------

SPECIMEN COLLECTION AND PREPARATION

Either citrate or EDTA collected platelet poor plasma may be used for this assay. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines-Fifth Edition", CLSI Document H21-A5, Vol. 28, No. 5, 2008¹⁵. Plasma collection should be performed as follows:

- Collect 9 parts of blood into 1 part of 3.2% (0.109M) trisodium citrate anticoagulant solution.
- Centrifuge the blood sample at 1,500 x g for 15 minutes.
- Plasma should be stored at 2° to 8°C and assayed within 2 hours. Alternatively, plasma may be stored for:

	2 weeks	6 months
	-20°C	-70°C

- Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at 2° to 8°C and assayed within 2 hours.

ASSAY PROCEDURE

- Six (6) antibody coated microwell strips are provided assembled in a plastic holder, giving the appearance of a fixed 96-microwell plate. Please note that the strips are constructed such that they fit into the frame in only one orientation. Open the foil pouch and remove the strip/frame assembly. In the foil pouch is a packet of desiccant. Do not discard this desiccant. Keep the desiccant in the pouch for storing unused microwell strips. Determine the number of microwells needed to generate the standard curve and assay the number of unknown specimen samples and references/ controls. Remember, it is recommended to assay standards and samples in duplicate. Remove the strips that will not be used, replace in the foil pouch and tightly reseal the pouch with the desiccant inside. Store the sealed pouch at 2° - 8°C until the expiration date printed on the pouch label.
- Add 50 µL of prepared PET Buffer to each microwell.
- Dilution of plasma samples is seldom necessary, but samples containing more than 50 ng/mL PAI-1 can be diluted with the 0 ng/mL PAI-1 Standard or tPA/PAI Depleted Plasma (BioMedica Diagnostics REF 273). The amounts of PAI-1 are then deduced from the calibration curve by multiplying the measured sample value by the dilution factor.
- Add 20 µL of PAI-1 standard or plasma sample to each microwell. Cover the wells with the acetate sheet provided. Place the strip/frame assembly on an orbital microwell plate shaker with agitation (at 500-600 rpm) and incubate for 1 hour at room temperature (20°-25°C).
- Remove the acetate cover sheet (do not discard) and add 50 µL of the Detection Antibody to each microwell. **Do not empty the wells.** Replace the acetate sheet over the microwells and incubate on the plate shaker (at 500-600 rpm) for 1 hour at room temperature (20°-25°C).
- Remove the acetate sheet and empty the contents of the microwells. Wash the wells four times with the PET Buffer. Washing can either be performed using microwell plate washing equipment or manually (fill the wells with PET Buffer with a pipette or squeeze bottle, wait three minutes, empty and remove droplets by tapping the plate 4-5 times face down against absorbing material. Repeat four times.
- Add Substrate as follows. If using all 96 microwells: add 30 µL of 30% H₂O₂ to all 24 mL of Substrate immediately before use and mix well. If using fewer than all 96 microwells: for each 16 microwells used, add 5 µL of 30% H₂O₂ to a 4.0 mL aliquot of Substrate and mix well. The H₂O₂ should be added to the Substrate just prior to use. Add 200 µL of the Substrate/H₂O₂ mixture to each microwell. Replace the acetate sheet over the microwells and incubate on the plate shaker (at 500-600 rpm) for 12 minutes at room temperature (20°-25°C).
- Stop the enzymatic reaction by adding 50 µL of Stop Solution. Add Stop Solution with the same speed and order as you added the Substrate. Store the strips in the dark until measuring the absorbances.
- Read the absorbance of the Substrate solutions at 490 nm on a microwell plate reader within 15 minutes.


Procedural Note

Completely fill the microwells during the wash step and make sure the microwells are completely empty after each wash. Do not allow the microwells to dry out.

RESULTS

A standard curve is constructed by plotting the mean absorbance value measured for each PAI-1 standard versus its corresponding PAI-1 concentration. A standard curve should be constructed each time the assay is performed. The following curve is presented for demonstration purposes only.

Representative Standard Curve

CALCULATIONS

Interpolate the PAI-1 concentration of the plasma sample directly from the standard curve. A curve regression software function that is included with the microwell plate reader may be used to calculate the concentrations.

Plasma samples found to contain more than 50 ng/mL PAI-1 should be diluted with the 0 ng/mL PAI-1 Standard supplied or tPA/PAI Depleted Plasma (BioMedica Diagnostics REF 273).

TRACEABILITY OF CONTROL MATERIALS

Information regarding traceability of control materials is available upon request from BioMedica Diagnostics.¹⁶

QUALITY CONTROL

It is recommended that 2 plasma samples containing PAI-1 antigen between 6 to 14 ng/mL (low control) and 36 to 44 ng/mL (high control), be stored at -20°C or below, in small aliquots and used as quality control standards each time the assay is run. Failure to obtain a PAI-1 antigen level within two standard deviations of the mean for each control standard may invalidate the assay. The 50 ng/mL standard and the 0 ng/mL standard are provided to construct the standard curve and must not be used to check device performance.

LIMITATIONS OF PROCEDURE

Many studies have shown that platelets are a source of PAI-1 in the bloodstream.^{1,2} Specimens not prepared according to the instructions, may be contaminated with platelets and therefore may have artificially high PAI-1 values. DO NOT subject plasma samples to more than one freeze-thaw cycle. As centrifugation cannot generate platelet-depleted plasma, there will be some platelets present. Successive freeze-thaw cycles will lyse any platelets and release PAI-1 into the sample. Plasmas subjected to more than one freeze-thaw cycle may have artificially high PAI-1 measurements.

Sodium azide present in a test sample will inhibit the peroxidase activity of the detection antibody. If azide was present in the sample, wash the well after the sample incubation and refill the well with 70 µL PET buffer before adding the conjugate.

Bacterial contamination in the filtered deionised water may cause aberrant results.

Each laboratory must determine separate normal ranges for EDTA or citrated collected plasma samples.

INTERFERENCES

The presence of bilirubin (1.5 mg/dL), cholesterol (228 mg/dL) and the appearance of slight to moderate hemolysis do not influence the determination of human PAI-1 in plasma.

EXPECTED VALUES

In a group of 167 healthy adults, 95% of plasma PAI-1 values were between 2 and 47 ng/mL using citrated plasmas. Each laboratory must determine its own normal range for both citrated and EDTA collected plasmas.

Studies show that PAI-1 activity correlates well with PAI-1 antigen in plasma samples.^{1,2} Nonetheless, platelets contain stored PAI-1 that is released in an inactive form, affecting PAI-1 antigen values but not PAI-1 activity values, so platelet release should be avoided. Values of PAI-1 antigen in human platelet-poor plasma have been reported to range between 4 and 43 ng/mL (mean of 18 ± 10 ng/mL with a correlation to PAI-1 activity of $r=0.80$).² Patients with recurrent deep vein thrombosis had higher PAI-1 antigen levels (mean of 44 ± 20 ng/mL).² PAI-1 antigen levels are also elevated during the third trimester of normal pregnancy.^{2, 11}

PERFORMANCE CHARACTERISTICS

Precision

Studies evaluated the intra-assay and inter-assay variations of two sample concentrations assayed in quadruplicate over 10 runs (n=40). The following variations were found.¹⁶

Sample	Mean Assayed PAI-1 Value	Intra-Assay Variation	Inter-Assay Variation
No. 1	19 ng/mL	6.6%	9.0%
No. 2	40 ng/mL	5.4%	6.9%

Sensitivity

The assay quantifies various forms of human Plasminogen Activator Inhibitor Type-1 (PAI-1) antigen. It detects active and inactive forms of PAI-1 as well as PAI-1 complexed with tissue-type Plasminogen Activator (tPA) and urokinase-type Plasminogen Activator (uPA). The lower limit detection is 2.2 ng/mL PAI-1 in the undiluted sample when the assay is performed according to the procedure.¹⁶

REFERENCES

1. Sprengers, E.D., *et al.* Plasminogen activator inhibitor. *Blood* 1987; **69**: 381-87.
2. Declerck, P.J., *et al.* Measurement of plasminogen activator inhibitor 1 (PAI-1) in biological fluids with a murine monoclonal antibody based enzyme-linked immunosorbent assay. *Blood* 1988; **71**: 220-225.
3. Schleef, R.R., *et al.* Bleeding diathesis due to decreased functional activity of type 1 plasminogen activator inhibitor. *Journal of Clinical Investigation* 1989; **83**: 1747.
4. Dieval, J., *et al.* A lifelong bleeding disorder associated with a deficiency of plasminogen activator inhibitor type 1. *Blood* 1991; **77**: 528.
5. Hamsten, A., *et al.* Plasminogen activator inhibitor in plasma: Risk factor for recurrent myocardial infarction. *Lancet* 1987; **2**: 3-8.
6. Juhan-Vague, I., *et al.* On Behalf of the ECAT Study Group. Fibrinolytic Factors and the Risk of Myocardial Infarction or Sudden Death in Patients with Angina Pectoris. *Circulation* 1996; **94**: 2057-2063.
7. Thøgersen, A.M., *et al.* High plasminogen activator inhibitor and tissue plasminogen activator levels in plasma precede a first acute myocardial infarction in both men and women: evidence for the fibrinolytic system as an independent primary risk factor. *Circulation* 1998; **98**: 2241-47.
8. Bergstein, J.M., *et al.* Role of plasminogen activator inhibitor type-1 in the pathogenesis and outcome of the Hemolytic Uremic Syndrome. *New England Journal of Medicine* 1992; **327**: 755-759.
9. Margaglione, M., *et al.* Abnormally High Circulation Levels of Tissue Plasminogen Activator and Plasminogen Inhibitor-1 in Patients with a History of Ischemic Stroke. *Arteriosclerosis and Thrombosis* 1994; **14**: 1741-1745.
10. He, S., *et al.* Increased blood flow resistance in placental circulation and levels of plasminogen activator inhibitors type 1 and 2 in severe preeclampsia. *Blood, Coagulation & Fibrinolysis* 1995; **6**: 703-708.
11. Cerneca, F., *et al.* Coagulation and fibrinolysis changes in normal pregnancy: Increased levels of procoagulants and reduced levels of inhibitors during pregnancy include a hypercoagulable state, combined with a reactive fibrinolysis. *Obstetrics & Gynecology* 1997; **73**: 31-36.
12. Bellart, J., *et al.* Coagulation and fibrinolysis parameters in normal pregnancy and in gestational diabetes. *American Journal of Perinatology* 1998; **15**: 479-486.
13. Tran-Thang, C., *et al.* Plasminogen activators and plasminogen activator inhibitors in liver deficiencies caused by chronic alcoholism or infectious hepatitis. *Thrombosis & Haemostasis* 1989; **62**: 651.
14. EC Regulation No. 1272/2008.
15. Collection, Transport and Processing of Blood Specimens for Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. *CLSI Document H21-A5*, 2008; Vol. 28, No. 5.
16. Data on file at BioMedica Diagnostics, Stamford, Connecticut, USA 06902.

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LITERATURHINWEISE

1. Sprengers, E.D., et al. Plasminogen activator inhibitor. *Blood* 1987; **69**: 381-87.
2. Declerck, P.J., et al. Measurement of plasminogen activator inhibitor 1 (PAI-1) in biological fluids with a murine monoclonal antibody based enzyme-linked immunosorbent assay. *Blood* 1988; **71**: 220-225.
3. Sleaf, R.R., et al. Bleeding diathesis due to decreased functional activity of type 1 plasminogen activator inhibitor. *Journal of Clinical Investigation* 1989; **83**: 1747.

4. Dieval, J., et al. A lifelong bleeding disorder associated with a deficiency of plasminogen activator inhibitor type 1. *Blood* 1991; **77**: 528.
5. Hamsten, A., et al. Plasminogen activator inhibitor in plasma: Risk factor for recurrent myocardial infarction. *Lancet* 1987; **2**: 3-8.
6. Juhan-Vague, I., et al. on behalf of the ECAT Study Group. Fibrinolytic Factors and the Risk of Myocardial Infarction or Sudden Death in Patients with Angina Pectoris. *Circulation* 1996; **94**: 2057-2063.
7. Thøgersen, A.M., et al. High plasminogen activator inhibitor and tissue plasminogen activator levels in plasma precede a first acute myocardial infarction in both men and women: evidence for the fibrinolytic system as an independent primary risk factor. *Circulation* 1998; **98**: 2241-47.
8. Bergstein, J.M., et al. Role of plasminogen activator inhibitor type-1 in the pathogenesis and outcome of the Hemolytic Uremic Syndrome. *New England Journal of Medicine* 1992; **327**, 755-759.
9. Margaglione, M., et al. Abnormally High Circulation Levels of Tissue Plasminogen Activator and Plasminogen Inhibitor-1 in Patients with a History of Ischemic Stroke. *Arteriosclerosis and Thrombosis* 1994; **14**: 1741-1745.
10. He, S., et al. Increased blood flow resistance in placental circulation and levels of plasminogen activator inhibitors type 1 and 2 in severe preeclampsia. *Blood, Coagulation & Fibrinolysis* 1995; **6**: 703-708.
11. Cerneca, F., et al. Coagulation and fibrinolysis changes in normal pregnancy: Increased levels of procoagulants and reduced levels of inhibitors during pregnancy include a hypercoagulable state, combined with a reactive fibrinolysis. *Obstetrics & Gynecology* 1997; **73**: 31-36.
12. Bellart, J., et al. Coagulation and fibrinolysis parameters in normal pregnancy and in gestational diabetes. *American Journal of Perinatology* 1998; **15**: 479-486.
13. Tran-Thang, C., et al. Plasminogen activators and plasminogen activator inhibitors in liver deficiencies caused by chronic alcoholism or infectious hepatitis. *Thrombosis & Haemostasis* 1989; **62**: 651.
14. EG Richtlinie 1272/2008.
15. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. *CLSI Document H21-A5*, 2008; Vol. 28, No. 5.
16. Daten liegen bei BioMedica Diagnostics, Stamford, Connecticut, USA 06902.

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REFERENCES

1. Sprengers, E.D., *et al.* Plasminogen activator inhibitor. *Blood* 1987; **69**: 381-87.
2. Declerck, P.J., *et al.* Measurement of plasminogen activator inhibitor 1 (PAI-1) in biological fluids with a murine monoclonal antibody based enzyme-linked immunosorbent assay. *Blood* 1988; **71**: 220-225.
3. Schleef, R.R., *et al.* Bleeding diathesis due to decreased functional activity of type 1 plasminogen activator inhibitor. *Journal of Clinical Investigation* 1989; **83**: 1747.
4. Dieval, J., *et al.* A lifelong bleeding disorder associated with a deficiency of plasminogen activator inhibitor type 1. *Blood* 1991; **77**: 528.
5. Hamsten, A., *et al.* Plasminogen activator inhibitor in plasma: Risk factor for recurrent myocardial infarction. *Lancet* 1987; **2**: 3-8.
6. Juhan-Vague, I., *et al.* On Behalf of the ECAT Study Group. Fibrinolytic Factors and the Risk of Myocardial Infarction or Sudden Death in Patients with Angina Pectoris. *Circulation* 1996; **94**: 2057-2063.
7. Thögersen, A.M., *et al.* High plasminogen activator inhibitor and tissue plasminogen activator levels in plasma precede a first acute myocardial infarction in both men and women: evidence for the fibrinolytic system as an independent primary risk factor. *Circulation* 1998; **98**: 2241-47.
8. Bergstein, J.M., *et al.* Role of plasminogen activator inhibitor type-1 in the pathogenesis and outcome of the Hemolytic Uremic Syndrome. *New England Journal of Medicine* 1992; **327**: 755-759.
9. Margaglione, M., *et al.* Abnormally High Circulation Levels of Tissue Plasminogen Activator and Plasminogen Inhibitor-1 in Patients with a History of Ischemic Stroke. *Arteriosclerosis and Thrombosis* 1994; **14**: 1741-1745.
10. He, S., *et al.* Increased blood flow resistance in placental circulation and levels of plasminogen activator inhibitors type 1 and 2 in severe preeclampsia. *Blood, Coagulation & Fibrinolysis* 1995; **6**: 703-708.
11. Cerneca, F., *et al.* Coagulation and fibrinolysis changes in normal pregnancy: Increased levels of procoagulants and reduced levels of inhibitors during pregnancy include a hypercoagulable state, combined with a reactive fibrinolysis. *Obstetrics & Gynecology* 1997; **73**: 31-36.
12. Bellart, J., *et al.* Coagulation and fibrinolysis parameters in normal pregnancy and in gestational diabetes. *American Journal of Perinatology* 1998; **15**: 479-486.
13. Tran-Thang, C., *et al.* Plasminogen activators and plasminogen activator inhibitors in liver deficiencies caused by chronic alcoholism or infectious hepatitis. *Thrombosis & Haemostasis* 1989; **62**: 651.
14. EC Directive 1272/2008.
15. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. *CLSI Document H21-A5*, 2008; Vol. 28, No. 5.
16. Données sur le fichier, BioMedica Diagnostics, Stamford, Connecticut USA 06902.

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1. Sprengers, E.D., *et al.* Plasminogen activator inhibitor. *Blood* 1987; **69**: 381-87.
2. Declerck, P.J., *et al.* Measurement of plasminogen activator inhibitor 1 (PAI-1) in biological fluids with a murine monoclonal antibody based enzyme-linked immunosorbent assay. *Blood* 1988; **71**: 220-225.
3. Schleef, R.R., *et al.* Bleeding diathesis due to decreased functional activity of type 1 plasminogen activator inhibitor. *Journal of Clinical Investigation* 1989; **83**: 1747.
4. Dieval, J., *et al.* A lifelong bleeding disorder associated with a deficiency of plasminogen activator inhibitor type 1. *Blood* 1991; **77**: 528.
5. Hamsten, A., *et al.* Plasminogen activator inhibitor in plasma: Risk factor for recurrent myocardial infarction. *Lancet* 1987; **2**: 3-8.
6. Juhan-Vague, I., *et al.* on behalf of the ECAT Study Group. Fibrinolytic Factors and the Risk of Myocardial Infarction or Sudden Death in Patients with Angina Pectoris. *Circulation* 1996; **94**: 2057-2063.
7. Thøgersen, A.M., *et al.* High plasminogen activator inhibitor and tissue plasminogen activator levels in plasma precede a first acute myocardial infarction in both men and women: evidence for the fibrinolytic system as an independent primary risk factor. *Circulation* 1998; **98**: 2241-47.
8. Bergstein, J.M., *et al.* Role of plasminogen activator inhibitor type-1 in the pathogenesis and outcome of the Hemolytic Uremic Syndrome. *New England Journal of Medicine* 1992; **327**, 755-759.
9. Margaglione, M., *et al.* Abnormally High Circulation Levels of Tissue Plasminogen Activator and Plasminogen Inhibitor-1 in Patients with a History of Ischemic Stroke. *Arteriosclerosis and Thrombosis* 1994; **14**: 1741-1745.
10. He, S., *et al.* Increased blood flow resistance in placental circulation and levels of plasminogen activator inhibitors type 1 and 2 in severe preeclampsia. *Blood, Coagulation & Fibrinolysis* 1995; **6**: 703-708.
11. Cerneca, F., *et al.* Coagulation and fibrinolysis changes in normal pregnancy: Increased levels of procoagulants and reduced levels of inhibitors during pregnancy include a hypercoagulable state, combined with a reactive fibrinolysis. *Obstetrics & Gynecology* 1997; **73**: 31-36.
12. Bellart, J., *et al.* Coagulation and fibrinolysis parameters in normal pregnancy and in gestational diabetes. *American Journal of Perinatology* 1998; **15**: 479-486.
13. Tran-Thang, C., *et al.* Plasminogen activators and plasminogen activator inhibitors in liver deficiencies caused by chronic alcoholism or infectious hepatitis. *Thrombosis & Haemostasis* 1989; **62**: 651.
14. Directiva 1272/2008/EC
15. Collection, Transport and Preparation of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition. *CLSI Document H21-A5*, 2008; Vol. 28, No. 5.
16. Datos registrados a disposición de BioMedica Diagnostics, Stamford, Connecticut, USA 06902.

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BIBLIOGRAFIA

1. Sprengers, E.D., *et al.* Plasminogen activator inhibitor. *Blood* 1987; **69**: 381-87.
2. Declerck, P.J., *et al.* Measurement of plasminogen activator inhibitor 1 (PAI-1) in biological fluids with a murine monoclonal antibody based enzyme-linked immunosorbent assay. *Blood* 1988; **71**: 220-225.
3. Schleef, R.R., *et al.* Bleeding diathesis due to decreased functional activity of type 1 plasminogen activator inhibitor. *Journal of Clinical Investigation* 1989; **83**: 1747.

4. Dieval, J., *et al.* A lifelong bleeding disorder associated with a deficiency of plasminogen activator inhibitor type 1. *Blood* 1991; **77**: 528.
5. Hamsten, A., *et al.* Plasminogen activator inhibitor in plasma: Risk factor for recurrent myocardial infarction. *Lancet* 1987; **2**: 3-8.
6. Juhan-Vague, I., *et al.* on behalf of the ECAT Study Group. Fibrinolytic Factors and the Risk of Myocardial Infarction or Sudden Death in Patients with Angina Pectoris. *Circulation* 1996; **94**: 2057-2063.
7. Thøgersen, A.M., *et al.* High plasminogen activator inhibitor and tissue plasminogen activator levels in plasma precede a first acute myocardial infarction in both men and women: evidence for the fibrinolytic system as an independent primary risk factor. *Circulation* 1998; **98**: 2241-47.
8. Bergstein, J.M., *et al.* Role of plasminogen activator inhibitor type-1 in the pathogenesis and outcome of the Hemolytic Uremic Syndrome. *New England Journal of Medicine* 1992; **327**, 755-759.
9. Margaglione, M., *et al.* Abnormally High Circulation Levels of Tissue Plasminogen Activator and Plasminogen Inhibitor-1 in Patients with a History of Ischemic Stroke. *Arteriosclerosis and Thrombosis* 1994; **14**: 1741-1745.
10. He, S., *et al.* Increased blood flow resistance in placental circulation and levels of plasminogen activator inhibitors type 1 and 2 in severe preeclampsia. *Blood, Coagulation & Fibrinolysis* 1995; **6**: 703-708.
11. Cerneca, F., *et al.* Coagulation and fibrinolysis changes in normal pregnancy: Increased levels of procoagulants and reduced levels of inhibitors during pregnancy include a hypercoagulable state, combined with a reactive fibrinolysis. *Obstetrics & Gynecology* 1997; **73**: 31-36.
12. Bellart, J., *et al.* Coagulation and fibrinolysis parameters in normal pregnancy and in gestational diabetes. *American Journal of Perinatology* 1998; **15**: 479-486.
13. Tran-Thang, C., *et al.* Plasminogen activators and plasminogen activator inhibitors in liver deficiencies caused by chronic alcoholism or infectious hepatitis. *Thrombosis & Haemostasis* 1989; **62**: 651.
14. Direttiva CE 1272/2008.
15. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and of Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. *CLSI Document H21-A5*, 2008; Vol. 28, No. 5.
16. Dati su file disponibili presso la BioMedica Diagnostics, Stamford, Connecticut, USA 06902.

	ENGLISH	DEUTSCH	FRANCAIS	ESPAÑOL	ITALIANO
	Definition of Symbols	Bedeutung der Symbole	Definizione dei simboli	Definición de símbolos	Définition des symboles
	Consult instructions for use	Gebrauchsanleitung beachten	Lire le mode d'emploi	Consulte las instrucciones para el Uso	Consultare le istruzioni per l'uso
	Refer to SDS	SDS konsultieren	Reportez-vous à SDS	Consulte SDS	Consultare la SDS
	In vitro diagnostic medical device	Medizinisches Produkt zur in vitro-Diagnostik	Dispositif de diagnostic médical in vitro	Dispositivo para diagnóstico in vitro	Dispositivo medico per diagnosi in vitro
	Manufactured by	Hergestellt durch	Fabricant	Fabricado por	Prodotto da
	Temperature limitation: Store at 2°C to 8°C	Bei 2°C bis 8°C lagern	Conserver à une température de 2°C à 8°C	Conservar entre 2°C y 8°C	Conservare tra 2°C e 8°C
	Batch code / Lot number	Charge	Numéro de lot	Código de lote	Numero di lotto
	Expiration Date	Verfallsdatum	Date d'expiration	Fecha de caducidad	Data di scadenza
	Catalog number	Katalog-Nr	Numéro de référence dans le catalogue	Número de catálogo	Numero di catalogo
	CE Mark	CE-Siegel	Marque CE	Marca CE	Marchio CE
	European Authorised Representative	Autorisierte Vertretung für Europa	Représentant européen autorisé	Representante Europeo Autorizado	Rappresentante autorizzato per l'Europa
	Contains sufficient for <n> tests	Inhalt ausreichend für <n> Tests	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" tests
	Contains...	Enthält...	Contient...	Contiene...	Contiene...