
IMUBIND® ADAMTS13 Autoantibody ELISA

REF 814

EC

REP

Obelis s.a
Bd. Général Wahis 53, 1030 Brussels, BELGIUM

INTENDED USE

The IMUBIND® ADAMTS13 Autoantibody ELISA is intended for the measurement of ADAMTS13 IgG autoantibodies in human plasma. The assay is intended for in vitro diagnostic use.

EXPLANATION OF THE TEST

ADAMTS13, also known as von Willebrand Factor (vWF) Cleaving Protease, is a zinc metalloproteinase that cleaves ultra large vWF multimers (UL-vWF) at the Tyr(1605)-Met(1606) bond located in the A2 region of vWF.¹ Studies have shown that low levels of ADAMTS13 activity are associated with Thrombotic Thrombocytopenia Purpura (TTP), a life-threatening hematological condition characterized by low platelet count, microvascular thrombi, red cell fragmentation, CNS and renal complications.^{2,3} A deficiency or low level of ADAMTS13 activity (<5%) may lead to an accumulation of UL-vWF multimers.⁴ The UL-vWF multimers will bind to receptors on platelets inducing platelet aggregation and formation of intravascular thrombi.

Congenital TTP is a rare inherited disorder resulting from mutations in the ADAMTS13 gene. Genetic alterations have been identified at many sites within the ADAMTS13 gene and result in the production of non-functional ADAMTS13 protein.^{5,6} The acquired form of TTP is an autoimmune-like disorder caused by the development of autoantibodies to ADAMTS13 that inhibit its activity.⁷ The development of ADAMTS13 autoantibodies has been associated with certain drugs such as ticlopidine, clopidogrel, quinine and cyclosporine.

PRINCIPLE OF THE METHOD

Diluted plasma samples are added to microwells coated with a full-length recombinant ADAMTS13. During an incubation period, autoantibodies to ADAMTS13 present in the sample will bind to the protein coated to the wells. Following a washing step, a goat anti-human IgG antibody labeled with horseradish peroxidase (HRP) is added to the microwells and binds to IgG type ADAMTS13 autoantibodies that have bound to the coated microwells. Following another washing step, the addition of a perborate-3,3'-5,5'-tetramethylbenzidine (TMB) substrate and its subsequent reaction with the HRP present generates a blue colored solution. The reaction is stopped by adding sulfuric acid, which turns the solution color yellow. Measuring the solution absorbance at 450 nm and extrapolating the value with those of a standard curve determines the level of ADAMTS13 autoantibodies in the diluted plasma sample. ADAMTS13 autoantibody levels are expressed in Arbitrary Units/mL (AU/mL). One AU/mL is equal to 1 µg/mL of affinity purified human anti-ADAMTS13 IgG.

REAGENTS

- R1 ADAMTS13 coated microwells**, 96, with acetate cover sheet
R2 Assay Buffer, 15 mL, 2 vials (lyophilized)
R3 Plasma Standard, 60 AU/mL, 0.45 mL, 2 vials, (lyophilized)
R4 Positive Control, 0.35 mL, 2 vials (lyophilized)
R5 Detection Antibody, goat anti-human IgG-HRP, 140 µL
R6 Substrate, TMB, 11 mL
R7 Wash Buffer, PBS w/0.05% Tween 20, pH 7.4, 1 Liter (lyophilized)

WARNING

Source material for some of the reagents in this kit is of human origin. This material has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using FDA approved methods. As no known test method provides complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, reagents should be handled as recommended for any potentially infectious human specimen. Discard all waste associated with test specimens and human source reagents in a biohazard waste container.

Limited for research use only in the United States. Not for internal use in humans or animals. Do not use the kit components beyond the stated expiration date. Do not mix reagents from different kits. Avoid microbial contamination of the reagents. Do not smoke, eat or drink in areas in which specimens or kit reagents are handled. Do not pipette reagents by mouth. Wear laboratory coat and disposable gloves throughout the test procedure and wash hands thoroughly afterwards. Avoid splashing or aerosol formation.

Assay Buffer	Danger (US) Warning (EU)	 	CONT Polyethylene glycol octylphenol ether H317, H319, H334, P261, P264, P280, P302 + P352, P305 + P351 + P338, P333 + P313, P337 + P313
Detection Antibody	Warning		H319, P264, P280, P305 + P351 + P338, P337 + P313
Wash buffer, PBS + Tween	Warning		H319, P264, P280, P305 + P351 + P338, P337 + P313

Hazard Statements:	H317 May cause an allergic skin reaction. H319 Causes serious eye irritation. H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
Precautionary statements:	P261 Avoid breathing dust/fume/gas/mist/vapors/spray. P264 Wash thoroughly after handling. P280 Wear eye protection/face protection. Wear protective gloves. P302 + P352 If on skin: Wash with plenty of water. P305 + P351 + P338 If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P333 + P313 If skin irritation or rash occurs: Get medical advice/attention. P337 + P313 If eye irritation persists: Get medical advice/attention.

REAGENT PREPARATION AND STORAGE

Unopened and lyophilized reagents are stable until the expiration date printed on the box when stored as instructed.

- R1 ADAMTS13 Coated Microwells:** Once removed from the foil pouch, the microwells must be used within 30 minutes. Unused microwells may be stored at 2°-8°C for 4 weeks when sealed in the original pouch with the desiccant present, protected from any moisture.
- R2 Assay Buffer:** Reconstitute with 15 mL of filtered deionized/distilled water. Gently mix the contents of the vial. Assay Buffer may be used for up to 4 weeks when stored at 2°-8°C.
- R3 Plasma Standard:** Reconstitute with 0.45 mL of Assay Buffer and let stand for 10 minutes. Mix prior to use. Standard may be aliquoted and stored at -20°C for up to 6 months.
- R4 Positive Control:** Reconstitute with 0.35 mL of Assay Buffer. Once reconstituted, the control may be stored for up to 4 weeks at -20°C. The control must be diluted 1:2 or 1:4 with Assay Buffer prior to assaying.
- R5 Detection Antibody:** Prepare working strength Detection Antibody by diluting 1:100 with Assay Buffer. For running all 96 microwells at one time, dilute 100 µL of Detection Antibody Conjugate up to 10 mL in Assay Buffer. If all 96 microwells are not to be used, dilute 20 µL of Detection Antibody up to 2 mL in Assay Buffer for each 16 microwell strip that will be used. Working strength Detection Antibody may be stored in the dark for 4 hours at room temperature or for 72 hours at 2°-8°C.
- R6 Substrate, TMB:** Supplied ready to use. It may be used until the expiration date stated on the vial when stored in the dark at 2°-8°C.
- R7 Wash Buffer:** Dissolve the contents of the packet with 1 Liter of filtered deionized/distilled water. After reconstitution, the Wash Buffer may be used for up to 4 weeks when stored at 2°-8°C.

SPECIMEN COLLECTION AND PREPARATION

Only citrate collected platelet poor plasma may be used for this assay. Do Not Use EDTA. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines-Fifth Edition", CLSI Document H21-A5, Vol. 28, No. 5, 2008. Plasma collection should be performed as follows:

1. Collect 9 parts of blood into 1 part of 3.2% (0.109 M) trisodium citrate anticoagulant solution.
2. Centrifuge the blood sample at 10,000 x g for 15 minutes.
3. Plasma should be stored at 2°-8°C and assayed within 4 hours. Alternatively, plasma may be stored at -20°C for up to 6 months.
4. Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at 2°-8°C and assayed within 4 hours.

PROCEDURE

Materials Provided – See Reagents

Material Required But Not Provided

filtered deionized/distilled water

50-300 μ L eight channel multi-pipette

0-200 μ L, 200-1000 μ L single pipettes

microwell plate reader for reading absorbance at 450 nm

microwell plate orbital shaker, microwell plate washer (optional)

0.5 M H_2SO_4 **Caution:** Handle Sulfuric acid with great care. Avoid any skin and eye contact. Wear protection glasses and gloves when handling.

Preparing Standards

1. Open the foil pouch and remove the microwell strips/frame assembly. Remove the strips that will not be used, return them to the foil pouch and tightly reseal the pouch with the desiccant inside. Store the foil pouch at 2°-8°C.
2. Reconstitute the Plasma Standard as instructed under REAGENT PREPARATION. Prepare five (5) serial dilutions of the standard using one of the following methods (Step 3a or 3b below). The concentrations of the serially diluted standard will be 60, 30, 15, 7.5, 3.75 and 1.87 AU/mL, respectively.
- 3a. **Well Dilution Method:** Add 200 μ L of Standard to microwells A1/A2. Add 100 μ L of Assay Buffer to microwells B1/B2 – F1/F2. Serially dilute the Standard by pipetting 100 μ L of the standard from microwells A1/A2 into microwells B1/B2. Mix and pipette 100 μ L from wells B1/B2 to wells C1/C2. Repeat this process through wells F1/F2. Remove and discard 100 μ L from wells F1/F2. Add 100 μ L of Assay Buffer to wells G1/G2 which serves as the 0 AU/mL standard.
- 3b. **Microcentrifuge Tube Dilution Method.** Pipette 225 μ L Assay Buffer into 5 appropriate sized tubes, which are labeled 30, 15, 7.5, 3.75 or 1.87 AU/mL. Add 225 μ L of Standard into the tube labeled 30 AU/mL and mix. Remove 225 μ L of the diluted Standard from the 30 AU/mL tube and add this to 15 AU/mL tube and mix. Continue this process until all the tubes have diluted Standard. Transfer 100 μ L of undiluted Standard to wells A1 and A2. Transfer 100 μ L of the diluted Standard from the 30 AU/mL tube to wells B1 and B2. Similarly, transfer the remainder of the tubes with diluted Standard to the successive wells of the plate. Add 100 μ L of Assay Buffer to wells G1 and G2, which serves as the 0 AU/mL standard.

Assay Procedure

4. Dilute each plasma sample 1:20 (1 part plasma + 19 parts Assay Buffer). Add 100 μ L of Positive Control or diluted sample to a microwell, cover with the acetate sheet and incubate for 2.5 hours at room temperature (18-25°C) on an orbital microwell plate shaker with agitation (at 500 rpm).
5. Empty the contents of the microwells and wash 4 times with Wash Buffer. Washing may be performed either using microwell plate washing equipment or manually (fill the wells with Wash Buffer with a pipette or squeeze bottle, wait three minutes, empty and remove droplets by tapping the plate 4-5 times face down against absorbing material).
6. Add 100 μ L of Detection Antibody to each microwell, cover with the acetate sheet and incubate the wells for 1 hour at room temperature (18-25°C) on an orbital microwell plate shaker with agitation (at 500 rpm).
7. Wash the wells by repeating Step 5.
8. Add 100 μ L of Substrate to each microwell immediately after the wash step, cover the wells with the acetate sheet and incubate for 10 minutes at room temperature (18°-25°C). A blue color will develop.
9. Stop the enzymatic reaction by adding 50 μ L of 0.5M H_2SO_4 to each microwell. Add the acid with the same speed and order as you added the substrate. Shake the microwells to ensure even distribution of the H_2SO_4 . The solution color will turn yellow. Read the absorbances on a microwell plate reader at a wavelength of 450 nm within 10 minutes.


QUALITY CONTROL

The Positive Control plasma should be included each time the assay is run. The Positive Control is diluted at either a 1:2 or 1:4 dilution just prior to assaying. Since control materials are used for the effective monitoring of the performance of a test, each laboratory should use these values to determine out-of-control situations using their own laboratory rules.

RESULTS

Construct a standard curve by plotting the mean absorbance value for each standard versus its corresponding concentration of ADAMTS13. A standard curve should be generated each time the assay is performed. The following standard curve is for demonstration purposes only.

Representative Standard Curve

CALCULATIONS

The level of ADAMTS13 IgG autoantibody in the plasma sample is interpolated directly from the standard curve. Although the sample has been diluted 1:20 in the assay procedure, do not adjust the interpolated value to account for this dilution.

If the value of the sample exceeds the highest standard, dilute the sample further with Assay Buffer and repeat the assay. Account for this further dilution when calculating the value of the neat sample.

LIMITATIONS OF THE PROCEDURE

Platelet contamination in plasma samples will interfere with the assay results. Plasma samples must be free of platelets in order to have a validated result. Exercise great care in minimizing disruption of the platelet pellet while recovering the platelet poor plasma. Samples should not be frozen and thawed more than two times.

Samples should not be collected with EDTA as the anticoagulant. Icteric, lipemic and hemolyzed samples may interfere with the assay.

EXPECTED VALUES

In a study of 119 normal donor plasmas (n=119), the mean level of ADAMTS13 IgG autoantibodies was found to be 9.0 AU/mL with a 95% reference limit of 26.3 AU/mL.⁸ Each laboratory should establish its own normal range using the local population.

PERFORMANCE CHARACTERISTICS

Precision

The intra-assay coefficient of variation (CV) for this ELISA has been found to be 6.7%.

Specificity

The detection antibody used is highly specific for human IgG.

TRACEABILITY OF CALIBRATORS AND CONTROL MATERIAL

Information on traceability of calibrators and control material is available upon request.

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