

## Instructions for use

# Cortisol Urine ELISA

**REF**

**MS E-5100**



**IVD**



## **Cortisol Urine ELISA**

### **INTENDED USE**

Competitive immunoenzymatic colorimetric method for quantitative determination of Free Cortisol concentration in Urine

### **CLINICAL SIGNIFICANCE**

Cortisol is a steroid hormone released from the adrenal cortex in response to an hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

Cortisol is bound, in plasma, from corticosteroid-binding globulin (CBG, transcortin), with high affinity, and from albumin. Only free cortisol is available to most receptors.

These normal endogenous functions are the basis for the physiological consequences of chronic stress - prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune / inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

The free cortisol fraction represents the metabolically active cortisol. In normal conditions, less than 1% it comes excrete in urines. In pathological conditions (syndrome of Cushing) the levels of free urinary cortisol are elevated, because the CBG don't bound the plasmatic cortisol in excess and it was removed with urines.

During pregnancy or estrogen-progestogen treatment an increase of plasmatic cortisol caused by an increment of the production of the transport protein, but the levels of free urinary cortisol results normal to indicate a correct surrenic functionality.

This test is very useful to estimate the real surrenic function, because is dose the free cortisol, it is the metabolically active form. Moreover the measurement of free urinary cortisol is the better parameter for the diagnosis of the Cushing's syndrome.

### **PRINCIPLE**

Urinary Cortisol (antigen) in the sample competes with horseradish-peroxidase Cortisol (enzyme-labelled-antigen) for binding onto the limited number of anti-Cortisol (antibody) sites on the microplates (solid phase). After incubation, the bound/free separation is performed by a simple solid-phase washing.

The enzyme substrate ( $H_2O_2$ ) and the TMB-Substrate (TMB) are added. After an appropriate time has elapsed for maximum colour development, the enzyme reaction is stopped and the absorbances are determined.

Urinary Cortisol concentration in the sample is calculated based on a series of standard.

The colour intensity is inversely proportional to the urinary cortisol concentration in the sample.

## **Reagent, material and instrumentation**

### **Reagent and material supplied in the kit**

#### **Standards**

|                   | <b>Cat. no.</b>  | <b>Standard</b> | <b>Concentration</b> | <b>Volume/Vial</b> |
|-------------------|------------------|-----------------|----------------------|--------------------|
| <b>STANDARD A</b> | <b>MS E-5101</b> | Standard 0      | 0 ng/ml              | 1 ml               |
| <b>STANDARD B</b> | <b>MS E-5102</b> | Standard 1      | 10 ng/ml             | 1 ml               |
| <b>STANDARD C</b> | <b>MS E-5103</b> | Standard 2      | 50 ng/ml             | 1 ml               |
| <b>STANDARD D</b> | <b>MS E-5104</b> | Standard 3      | 150 ng/ml            | 1 ml               |
| <b>STANDARD E</b> | <b>MS E-5105</b> | Standard 4      | 500 ng/ml            | 1 ml               |

SD0 – STD4, 5x (1 vial = 1 mL)

#### **INC-BUFF**

#### **MS E-5113 Incubation buffer**

(1 bottle) 100 mL; Phosphate buffer 50 mM pH 7.4; BSA 1 g/L

#### **CONJUGATE-CONC**

#### **MS E-5140 Conjugate**

(1 bottle) 1 mL; Cortisol-HRP conjugate

#### **96**

#### **MS E-5131 Coated Microplate**

(1 microplate breakable); Anti-Cortisol-IgG adsorbed on microplate

#### **SUBSTRATE**

#### **MS E-0055 TMB-substrate**

(1 bottle) 15 mL; H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (avoid any skin contact)

#### **STOP-SOLN**

#### **MS E-0080 Stop solution**

(1 bottle) 15 mL; Sulphuric acid 0.15 mol/L (avoid any skin contact)

#### **CONTROL 1**

#### **MS E-5151 Low Control**

(1 vial = 1 mL); Ready to use

#### **CONTROL 2**

#### **MS E-5152 High Control**

(1 vial = 1 mL); Ready to use

### **Reagents necessary not supplied**

Distilled water

### **Auxiliary materials and instrumentation**

Automatic dispenser

Microplates reader

### **Note**

Store all reagents at 2 °C – 8 °C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close immediately after use.

Do not remove the adhesive sheet from the unused strips.

## **PRECAUTION**

The reagent contains Proclin 300 as preservative.

Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants.

Maximum precision is required for reconstitution and dispensation of the reagents.

Do not use different lots of reagents.

Do not use heavily haemolized samples.

This method allows the determination of Cortisol from 10 ng/mL to 500 ng/mL.

The clinical significance of the Cortisol determination can be invalidated if the patient was treated with corticosteroids or natural or synthetic steroids.

## **PROCEDURE**

### **Preparation of the Standard (S<sub>0</sub>, S<sub>1</sub>, S<sub>2</sub>, S<sub>3</sub>, S<sub>4</sub>)**

The standards have the following concentration of Cortisol:

|       | S <sub>0</sub> | S <sub>1</sub> | S <sub>2</sub> | S <sub>3</sub> | S <sub>4</sub> |
|-------|----------------|----------------|----------------|----------------|----------------|
| ng/mL | 0              | 10             | 50             | 150            | 500            |

Stable when stored at +4 °C until the expiration date of the kit;  
once opened the standards are stable six months at +4 °C.

### **Preparation of diluted Conjugate**

Prepare immediately before use.

Add 10 µL Conjugate (reagent 3) to 1.0 mL of Incubation Buffer (reagent 2). Mix gently.

**Stable for 3 hours** at room temperature.

### **Preparation of the Sample**

The total volume of urine excreted during a 24 hours should be collected and mixed in a single container.

Urine samples which are not to be assayed immediately should be stored at 2 °C – 8 °C or at – 20°C for periods longer than a week.

**Dilute the urine sample (1 + 1) with Incubation buffer (e.g. 100 µL + 100 µL)**

### **Procedure**

As it is necessary to perform the determination in duplicate, prepare two wells for each of the four points of the standard curve (S<sub>0</sub>-S<sub>4</sub>), one for Blank.

| Reagent                                                                                                                                                                          | Standard/Control | Samples | Blank  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|---------|--------|
| Standard S <sub>0</sub> -S <sub>4</sub><br>Controls                                                                                                                              | 10 µL            |         |        |
| Diluted samples                                                                                                                                                                  |                  | 10 µL   |        |
| Diluted Conjugate                                                                                                                                                                | 300 µL           | 300 µL  |        |
| Incubate at 37°C for 1 hour.<br>Remove the contents from each well; wash the wells with 400 µL of distilled water. Repeat the washing procedure by draining the water completely |                  |         |        |
| TMB substrate                                                                                                                                                                    | 100 µL           | 100 µL  |        |
| Incubate at room temperature 22÷28°C for 15 minutes in the dark                                                                                                                  |                  |         |        |
| Stop solution                                                                                                                                                                    | 100 µL           | 100 µL  | 100 µL |
| Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank.                                                                                                    |                  |         |        |

## **QUALITY CONTROL**

Each laboratory should assay controls at normal, high and low levels range of Urinary Cortisol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## **LIMITATION OF PROCEDURE**

### **Assay Performance**

Sample(s), which are contaminated microbiologically, should not be used in the assay.

Highly lipemic or haemolysed specimen(s) should similarly not be used.

It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve.

Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time deviation during reaction.

Plate readers measure vertically. Do not touch the bottom of the wells.

Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.

### **Interpretation**

If computer controlled data reduction is used to calculate the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

## **RESULTS**

### **Mean Absorbance**

Calculate the mean of the absorbance ( $E_m$ ) for each point of the standard curve and of each sample

### **Standard Curve**

Plot the values of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points (e.g.: Four Parameter Logistic).

### **Calculation of Results**

As the standards are prediluted, the dilution factor 2 for the urine samples has not to be taken into account for the calculation. The concentration of the samples (in ng/mL) can be read directly from the standard curve. Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/mL.

To calculate the cortisol concentration in urine, calculate as above and correct for total volume of volume of urine collected in 24 hours:

$$\text{ng/mL} \times \text{Vol(mL) urine 24 h} / 1.000 = \mu\text{g Cortisol/24h}$$

### **Reference Value**

50 – 190  $\mu\text{g/24 hours}$

## **Performance and Characteristics**

### **Precision**

### **Intra-Assay Variation**

Within run variation was determined by replicate measurements (10x) of three different control sera in one assay. The within assay variability is  $\leq 7\%$ .

### **Inter-Assay Variation**

Between run variation was determined by replicate measurements (10x) of three different control sera in different lots of kit. The between assay variability is  $\leq 9\%$ .

### **Accuracy**

The recovery of 25 - 50 - 100 - 200 ng/mL of Urinary Cortisol added to a sample gave an average value ( $\pm$ SD) of  $101.89\% \pm 4.55\%$  with reference to the original concentrations.

### **Sensitivity**

The lowest detectable concentration of urinary cortisol that can be distinguished from the zero standard is 2.0 ng/ml at the 95 % confidence limit.

### **Specificity**

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

|                             |                        |
|-----------------------------|------------------------|
| Cortisol                    | 100%                   |
| Cortisone                   | 10.8%                  |
| 11 $\alpha$ -deoxycortisol  | 18.7%                  |
| Corticosterone              | 2.4%                   |
| Progesterone                | 0.1%                   |
| Aldosterone                 | $1 \times 10^{-2}\%$   |
| 11 $\alpha$ OH Progesterone | $1 \times 10^{-2}\%$   |
| Cholesterol                 | $< 1 \times 10^{-6}\%$ |

### **Correlation with RIA**

The Urinary Cortisol ELISA was compared to another commercially available Urinary Cortisol assay. 50 serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$y = 1.10x - 1.56$$

$$r^2 = 0.949$$

y = Urinary Cortisol ELISA

x = Cortisol Immunotech RIA

### **WASTE MANAGEMENT**

Reagents must be disposed off in accordance with local regulations.

### **BIBLIOGRAPHY**

1. Foster, L. B. and Dunn, R. T. Clin. Chem. 20/3 365(1974)
2. De Lacerda, et al J. Clin. Endocr. and Metab. 36,227 (1973)
3. Roller, E., et al Clin. Chim. Acta 66 319 (1976)
4. Kobayashi, Y., et al Steroids, 32 no 1(1978)
5. Akarawa, et al Anal. Biochem. 97 248 (1979)

## **TROUBLESHOOTING**

### **ERRORS / POSSIBLE CAUSES / SUGGESTIONS**

#### **No colorimetric reaction**

- no conjugate pipetted
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

#### **Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

#### **Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

#### **Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed)

#### **too high within-run CV%**

- reagents and/or strips not pre-warmed to room temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)

#### **too high between-run CV %**

- incubation conditions not constant (time, temperature)
  - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

---

### **Symbols:**

|                                                                                     |                              |                                                                                     |                  |                                                                                       |                                   |
|-------------------------------------------------------------------------------------|------------------------------|-------------------------------------------------------------------------------------|------------------|---------------------------------------------------------------------------------------|-----------------------------------|
|  | Storage temperature          |  | Manufacturer     |  | Contains sufficient for <n> tests |
|  | Expiry date                  |  | Batch code       |  | For in-vitro diagnostic use only! |
|  | Consult instructions for use |  | Content          |  | CE labelled                       |
|  | Caution                      |  | Catalogue number |  | For research use only!            |