Instructions for use **S 100B ELISA**

TM E-4800

S100B ELISA

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of S100B concentration in serum.

CLINICAL SIGNIFICANCE

S100 is a 20 kDa protein belonging to the S100/calmodulin/troponin C superfamily of EF-hand calcium-binding proteins. S100 was originally isolated from human brain and considered a glial-cell specific protein (1). Today, 20 monomers of the S100 family have been identified based on structural and functional similarities (2,3). Most of the S100 proteins exist as dimers and are expressed in a cell-specific manner. Two of the S100 monomers, designated S100 α and S100 β (4) are highly conserved between species and are found as homo-($\beta\beta$) and heterodimers ($\alpha\beta$) in central nervous system glial cells and in certain peripheral cells eg. Schwann cells, melanocytes, adipocytes, and chondrocytes (5). S100 $\alpha\beta$ and S100 $\beta\beta$ are also present in malignant tissues, most notably in melanoma and to a lesser extent in glioma, thyroid cell carcinoma and renal cell carcinoma (2). Determination of S100B (like S100 $\alpha\beta$ and S100 $\beta\beta$ units) in serum has been shown to be clinically useful for prognosis and treatment monitoring of patients diagnosed with malignant melanoma (6-9). Studies also suggest that S100B may be useful in the management of patients with brain damage from eg. traumatic head injury, perinatal asphyxia, cardiac arrest, cardiac surgery and stroke (10-13).

Principle

The S100B ELISA TEST is based on binding of S100B by two antibodies, one immobilized on microwell plates, and the other one conjugates with horseradish peroxidase (HRP).

The assay is on two steps binding procedure and after every incubation step, the bound/free separation is performed by a simple solid-phase washing, then the substrate solution (TMB) is added. After an appropriate time has elapsed for maximum color development, the enzyme reaction is stopped and the absorbances are determined. The S100B concentration in the sample is calculated based on a series of standard.

The color intensity is proportional to the S100B concentration in the sample.

Reagent, material and instrumentation

Reagent and material supplied in the kit.

Standards

	Cat. no.	Standard	Concentration	Volume/Vial
STANDARD A	TM E-4801	Standard 0	0 pg/mL	1 ml
STANDARD B	TM E-4802	Standard 1	80 pg/mL	1 ml
STANDARD C	TM E-4803	Standard 2	160 pg/mL	1 ml
STANDARD D	TM E-4804	Standard 3	320 pg/mL	1 ml
STANDARD E	TM E-4805	Standard 4	2000 pg/mL	1 ml
STANDARD F	TM E-4806	Standard 5	4500 pg/mL	1 ml

6x (1 vial = lyophilized); STD0 - STD5

CONTROL 1 + CONTROL 2 **TM E-4851 + TM E-4852 Controls**

2x (1 vial = lyophilized); Negative Control; Positive Control

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(1 vial) 0.6 mL; Anti-S100B-HRP conjugate

ASSAY-BUFF **TM E-4813 Assay buffer**

(1 bottle) 12 mL; Tris buffer; BSA 10 gr/L, stabilizing reagent

96 **TM E-4831 Coated microplate**

(1 microplate breakable); Anti-S100B adsorbed on microplate

SUBSTRATE **MS E-0055 TMB-Substrate**

(1 bottle) 12 mL; H_2O_2 -TMB 0.25 g/L (avoid any skin contact)

STOP-SOLN **MS E-0080 Stop solution**

(1 bottle) 12 mL; Sulphuric acid 0.15 mol/L (avoid any skin contact)

WASH-CONC 50x **SA E-0030 Conc.wash solution 50x**

(1 bottle) 20 mL; NaCl 45g/L Tween 20 55g/L

Necessary reagents not supplied with the kit

Distilled water.

Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader

Notes

Store all reagents between +2°C - 8°C in the dark.

Open the bag of reagent 6 (Coated Microplate) only when it is at room temperature and close immediately after use.

Do not remove the adhesive sheets on the strips unutilized

PRECAUTION

- The reagents contain Merthiolate Sodium as preservative.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- All the reagents have a lot to lot consistency; do not mix various lot numbers kits components within a test.
- This method allows the determination of S100B from 10 to 5000 pg/mL.
- The calibrator concentrations are lot specific and are reported on the vial labels.
- Do not use heavily hemolysed samples
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants.

PROCEDURE**Preparation of the sample**

The S100B determination must be carried out on serum. Do not use hemolyzed samples.

Samples can be stored at +2°C - 8°C for 1 day; for long periods store at -20°C. Avoid repeated freeze-thaw cycles. Do not allow the samples at room temperature for long period.

For sample with concentration higher than 5 ng/mL dilute the sample with Assay buffer

Preparation of the Standards and Controls

Reconstitute standards and controls with 1 mL of distilled water before use;

once reconstituted they are stable for 4 weeks at 2°C – 8°C and about six month if stored at -20°C.

It is advised to divide the content in aliquots and store them at -20°C.

The values of standard concentrations are reported on vial label.

Avoid repeated freeze-thaw cycles and long time exposure at room temperature.

Preparation of the Conjugate

Prepare 2 hour before use.

Add 50 µL Conjugate (reagent 4) to 1.0 mL of Conjugate Buffer (reagent 3).

(The quantity of diluted conjugate is proportional to the number of tests).

Mix gently for 5 minutes, with rotating mixer.

Stable for 3 hours at room temperature.

Preparation of the wash solution

Dilute 10 mL of Wash Solution Concentrate (50x) with 490 ml of distilled or deionized; for different volumes keep dilution ratio. Store at room temperature until the expiry date write on the wash solution concentrate label.

Procedure

As it is necessary to perform the determination in duplicate, prepare two wells for each point of the standard curve (S_0 - S_5), two for each Control, two for each sample, one for Blank.

Pipette:

Reagent	Standard	Sample	Blank
Standard S_0 - S_5 , Controls	50 μ L		
Sample		50 μ L	
Assay Buffer	50 μ L	50 μ L	
Incubate 2 h at room temperature ($22 \div 28^\circ\text{C}$). Remove the contents from each well, wash the wells three times with 300 μ L of diluted wash solution			
Diluted Conjugate	100 μ L	100 μ L	
Incubate 1 h at room temperature ($22^\circ\text{C} - 28^\circ\text{C}$). Remove the contents from each well, wash the wells three times with 300 μ L of diluted wash solution			
TMB-Substrate	100 μ L	100 μ L	100 μ L
Incubate 30 minutes at room temperature $22 \div 28^\circ\text{C}$, in the dark.			
Stop solution	100 μ L	100 μ L	100 μ L
Read the absorbance (E) at 450 nm against Blank			

RESULTS

Mean Absorbance

Calculate the mean of the absorbencies (E_m) for each point of the standard curve and of each sample.

Standard Curve

Plot the values of absorbance of the standards ($S_0 - S_5$) against concentration. Draw the best-fit curve through the plotted points (e.g.: Cubic spline or Four Parameter Logistic).

Calculation of Results

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations in pg/mL.

REFERENCES VALUES

Each laboratory must establish its own normal ranges based on patient population:

Normal range: Mean = 50 pg/mL; SD = 15 pg/mL

Pathological limit: > 75 pg/mL

Performance Characteristics

Specificity

The antibody recognizes specifically the β subunit therefore reacts with S100 $\alpha\beta$ and S100 $\beta\beta$ units.

It is not reactive against the S100 $\alpha\alpha$ unit.

Cross-reacts with S100 from bovine, porcine, rabbit cat and rat.

Does not react with the other EF-hand family proteins.

Sensitivity

The minimal detectable concentration of S100 β by this assay is estimated to be 10 pg/mL.

Limitations of Use:

In this assay, no hook effect is observed up to 5000 pg/mL of S100B

WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations

BIBLIOGRAPHY

1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins: structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T et al., (2003) Use of serum S-100 and S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E, et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicenter study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Stigbrand T, Nyberg L, Ullen A, Haglid K, Sandstrom E, Brundell J. A new specific method for measuring S-100B in serum. *Int J Biol Markers*. 2000 Jan-Mar;15(1):33-40.
15. Chen DQ, Zhu LL. Dynamic change of serum protein S100b and its clinical significance in patients with traumatic brain injury.. *Chin J Traumatol*, 2005 Aug 1;8(4):245-8.
16. Sawauchi S, Taya K, Murakami S, Ishi T, Ohtsuka T, Kato N, Kaku S, Tanaka T, Morooka S, Yuhki K, Urashima M, Abe T. Serum S-100B protein and neuron-specific enolase after traumatic brain injury. *No Shinkei Geka*. 2005 Nov;33(11):1073-80

TROUBLESHOOTING

ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers -insufficient washing (conjugates not properly removed)

too high within-run CV%

- reagents and/or strips not pre-warmed to Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)

too high between-run CV %

- incubation conditions not constant (time, temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled
	Caution		Catalogue number		For research use only!

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1. Moore BW (1965) A soluble protein characteristic of the nervous system. *Biochem Biophys Res Commun* 19:739-744.
2. Zimmer DB et al., (1995) The S100 protein family history, function and expression. *Brain Res Bull* 37:417-429.
3. Heizmann CW et al., (2002) S100 proteins: structure, functions and pathology. *Front Biosci* 7:1356-1368.
4. Schäfer BW et al. (1995) Isolation of a YAC clone covering a cluster of nine S100 genes on human chromosome 1q21: rationale for a new nomenclature of the S100 calcium-binding protein family. *Genomics* 25:638-643.
5. Takahashi K et al., (1984) Immunohistochemical study on the distribution of α and β subunits of S-100 protein in human neoplasm and normal tissues. *Virchows Arch* 45:385-396.
6. Banfalvi T et al., (2003) Use of serum S-100B protein levels to monitor the clinical course of malignant melanoma. *Eur J Cancer* 39:164-169.
7. Djureen-Mårtensson E, et al., (2001) Serum S-100b protein as a prognostic marker in malignant cutaneous melanoma. *J Clin Oncol* 19:824-831.
8. Hauschild A et al., (1999) S100B protein detection in serum is significant prognostic factor in metastatic melanoma. *Oncology* 56:338-344.
9. Wunderlich MT et al., (1999) Early Neurobehavioral Outcome after stroke is related to release of neurobiochemical markers of brain damage. *Stroke* 30:1190-1195.
10. Martens P et al., (1998) Serum S100 and neuron specific enolase for prediction of regaining consciousness after global cerebral ischemia. *Stroke* 29:2363-2366.
11. Rosén H et al., (1998) Increased serum levels of the S-100 protein are associated with hypoxic brain damage after cardiac arrest. *Stroke* 29: 473-477.
12. Ingebrigtsen T et al., (2000) The clinical value of serum S-100 protein measurements in minor head injury: a Scandinavian multicenter study. *Brain Inj* 14:1047-1055
13. Michetti F and Gazzolo D (2002) S100B protein in biological fluids: A tool for perinatal medicine *Clin Chem* 48:2097-2104
14. Stigbrand T, Nyberg L, Ullen A, Haglid K, Sandstrom E, Brundell J. A new specific method for measuring S-100B in serum. *Int J Biol Markers*. 2000 Jan-Mar;15(1):33-40.
15. Chen DQ, Zhu LL. Dynamic change of serum protein S100b and its clinical significance in patients with traumatic brain injury.. *Chin J Traumatol*, 2005 Aug 1;8(4):245-8.
16. Sawauchi S, Taya K, Murakami S, Ishi T, Ohtsuka T, Kato N, Kaku S, Tanaka T, Morooka S, Yuhki K, Urashima M, Abe T. Serum S-100B protein and neuron-specific enolase after traumatic brain injury. *No Shinkei Geka*. 2005 Nov;33(11):1073-80