

Borrelia ViraChip® IgM Test Kit

Chip-Immunoblot for the qualitative detection of **IgM** antibodies against specific **Borrelia species** antigens in human serum.

The **Borrelia ViraChip® IgM Test Kit** is an immunoblot based on an enzyme-immunoassay in a microarray format, carrying highly purified specific native antigens from *Borrelia afzelii* (Pko) and *Borrelia burgdorferi sensu stricto* as well as recombinant VlsE at defined positions. The Borrelia ViraChip® IgM Test Kit is manufactured according to the guidelines **98/79/EG** and **DIN 58967-40**.

The Borrelia ViraChip® IgM Test Kit fulfils the quality standards of the microbiological-infectiological guidelines ("**MiQ**" **12/2000**), the general DIN recommendation (**DIN 58967-40**) and can be evaluated as confirmatory assay according to "**MiQ**" **12/2000** (18) and **DIN 58969-44** (5) criteria. These guidelines postulate a *two-band* criterion for a positive result and describe special requirements for the detection of antibodies against *Borrelia burgdorferi*.

Principle of the assay

One microarray is fixed at the bottom of each well in a standard microtiter plate (MTP). The single breakable wells are stored in a holding frame with 96 positions. During the serum incubation step *Borrelia* specific IgM antibodies bind to the immobilised antigens, herein after referred to as spots, on the microarray. During the conjugate reaction, AP-conjugate binds to the antigen-antibody complex. The alkaline phosphatase converts the chromogen/substrate and thus, stains the antigen-antibody complex on the microarray purple. The washing procedures following serum, conjugate and chromogen/substrate incubation steps remove unbound antibodies and reagents.

The control spots include **serum controls, conjugate controls, calibrator controls** and a **negative control**. The analyte spots serve to detect antibodies against **p41, p39, OspC, Osp17** and **VlsE**.

Each well is coded by a colour system: The Borrelia ViraChip® IgM is marked with a **green quadrant**.

Order No.:	V-BSCMOK	Order No.:	V-BSCMDK (Deca Kit)
Kit size:	96 wells	Kit size:	10x 96 wells
Specimen:	10 µl serum	Specimen:	10 µl serum
Time for testing:	approx. 130 minutes	Time for testing:	approx. 130 minutes

Materials provided

1x or 10x 96 wells	Borrelia ViraChip® IgM Antigen Coated Wells	(Prod. No.: V-BSCMAC)
1x or 10x 1.5 ml	Wells with ViraChip® Microarrays, ready to use, green quadrant	
1x or 10x 100 ml	ViraChip® AP-Anti-Human IgM Conjugate	(Order No.: V-UVNMKI15)
1x or 10x 5 g	Universal anti-human IgM conjugate for ViraChip® tests	
1x or 10x 12 ml	ViraChip® / ViraStripe® / ViraBlot® Diluent / Wash Buffer	(Order No.: V-UVNUWP)
	Universal Diluent / Wash Buffer for ViraChip®, ViraStripe® and ViraBlot®	
	ViraChip® / ViraStripe® / ViraBlot® Diluent / Wash Powder	(Order No.: V-UVNUPM)
	Universal Diluent / Wash Powder for ViraChip®, ViraStripe® and ViraBlot®	
	ViraChip® Chromogen / Substrate Solution	(Order No.: V-UVUCUS)
	Universal Chromogen / Substrate Solution for ViraChip® tests	

Additionally available

330 µl	Borrelia ViraChip® IgM Positive Control	(Order No.: V-BSCMPK)
	Human, ready to use	
330 µl	Borrelia ViraChip® IgM,A,M Negative Control	(Order No.: V-BSCPNK)
	Human, ready to use	
96 wells	Microtiter plate with 96 empty wells	(Order No.: V-UVNMTP)

Preparation of reagents and patient samples

Bring all reagents to room temperature (20-23°C) prior to use. Let the packed microtiter plate acclimatise for at least 30 min before opening. The test has to be performed at room temperature and at a maximal humidity of 60%.

Diluent / Wash Buffer	Dilute Diluent / Wash Buffer Concentrate 1:10 with distilled or deionised water (100 ml concentrate + 900 ml water). Add Diluent / Wash Powder completely and stir well until all powder is dissolved. If needed, place onto a magnetic stirrer for 10-15 minutes. The pH value should be around pH 7.5 at 20°C.
Working Dilution:	
Wells:	Carefully unpack the microtiter plate (MTP) and place the required number of wells in an empty holding frame (see assay procedure, step 1). Use wells directly after removing from packing. Return unused wells directly into the original packing, seal accurately and store at 2-8°C.
Patient samples:	Use 100 µl of a 1:76 dilution of patient serum, e.g. 10 µl of patient serum plus 750 µl Diluent / Wash Buffer Working Dilution* .
Controls:	Optionally, use 100 µl of a 1:16 dilution of control serum, e.g. 10 µl of control serum plus 150 µl Diluent / Wash Buffer Working Dilution* .
Conjugate Working Dilution:	Prepare Conjugate Concentrate 1:10 with Diluent / Wash Buffer Working Dilution (see table 1). Prepare freshly prior to each test run. Do not store for further use.
Chromogen / Substrate Solution	Ready to use.

*) Depending on the equipment dilutions may be performed in several steps.

Borrelia ViraChip® IgM Test Kit

- 2 -

Preparation of Conjugate Working Dilution

Number of wells	Diluent / Wash Buffer Working Dilution		Conjugate Concentrate	Final Volume	Number of wells	Diluent / Wash Buffer Working Dilution		Conjugate Concentrate	Final Volume
1	0.090 ml	+	0.010 ml	0.10 ml	51	4.590 ml	+	0.510 ml	5.10 ml
2	0.180 ml	+	0.020 ml	0.20 ml	52	4.680 ml	+	0.520 ml	5.20 ml
3	0.270 ml	+	0.030 ml	0.30 ml	53	4.770 ml	+	0.530 ml	5.30 ml
4	0.360 ml	+	0.040 ml	0.40 ml	54	4.860 ml	+	0.540 ml	5.40 ml
5	0.450 ml	+	0.050 ml	0.50 ml	55	4.950 ml	+	0.550 ml	5.50 ml
6	0.540 ml	+	0.060 ml	0.60 ml	56	5.040 ml	+	0.560 ml	5.60 ml
7	0.630 ml	+	0.070 ml	0.70 ml	57	5.130 ml	+	0.570 ml	5.70 ml
8	0.720 ml	+	0.080 ml	0.80 ml	58	5.220 ml	+	0.580 ml	5.80 ml
9	0.810 ml	+	0.090 ml	0.90 ml	59	5.310 ml	+	0.590 ml	5.90 ml
10	0.900 ml	+	0.100 ml	1.00 ml	60	5.400 ml	+	0.600 ml	6.00 ml
11	0.990 ml	+	0.110 ml	1.10 ml	61	5.490 ml	+	0.610 ml	6.10 ml
12	1.080 ml	+	0.120 ml	1.20 ml	62	5.580 ml	+	0.620 ml	6.20 ml
13	1.170 ml	+	0.130 ml	1.30 ml	63	5.670 ml	+	0.630 ml	6.30 ml
14	1.260 ml	+	0.140 ml	1.40 ml	64	5.760 ml	+	0.640 ml	6.40 ml
15	1.350 ml	+	0.150 ml	1.50 ml	65	5.850 ml	+	0.650 ml	6.50 ml
16	1.440 ml	+	0.160 ml	1.60 ml	66	5.940 ml	+	0.660 ml	6.60 ml
17	1.530 ml	+	0.170 ml	1.70 ml	67	6.030 ml	+	0.670 ml	6.70 ml
18	1.620 ml	+	0.180 ml	1.80 ml	68	6.120 ml	+	0.680 ml	6.80 ml
19	1.710 ml	+	0.190 ml	1.90 ml	69	6.210 ml	+	0.690 ml	6.90 ml
20	1.800 ml	+	0.200 ml	2.00 ml	70	6.300 ml	+	0.700 ml	7.00 ml
21	1.890 ml	+	0.210 ml	2.10 ml	71	6.390 ml	+	0.710 ml	7.10 ml
22	1.980 ml	+	0.220 ml	2.20 ml	72	6.480 ml	+	0.720 ml	7.20 ml
23	2.070 ml	+	0.230 ml	2.30 ml	73	6.570 ml	+	0.730 ml	7.30 ml
24	2.160 ml	+	0.240 ml	2.40 ml	74	6.660 ml	+	0.740 ml	7.40 ml
25	2.250 ml	+	0.250 ml	2.50 ml	75	6.750 ml	+	0.750 ml	7.50 ml
26	2.340 ml	+	0.260 ml	2.60 ml	76	6.840 ml	+	0.760 ml	7.60 ml
27	2.430 ml	+	0.270 ml	2.70 ml	77	6.930 ml	+	0.770 ml	7.70 ml
28	2.520 ml	+	0.280 ml	2.80 ml	78	7.020 ml	+	0.780 ml	7.80 ml
29	2.610 ml	+	0.290 ml	2.90 ml	79	7.110 ml	+	0.790 ml	7.90 ml
30	2.700 ml	+	0.300 ml	3.00 ml	80	7.200 ml	+	0.800 ml	8.00 ml
31	2.790 ml	+	0.310 ml	3.10 ml	81	7.290 ml	+	0.810 ml	8.10 ml
32	2.880 ml	+	0.320 ml	3.20 ml	82	7.380 ml	+	0.820 ml	8.20 ml
33	2.970 ml	+	0.330 ml	3.30 ml	83	7.470 ml	+	0.830 ml	8.30 ml
34	3.060 ml	+	0.340 ml	3.40 ml	84	7.560 ml	+	0.840 ml	8.40 ml
35	3.150 ml	+	0.350 ml	3.50 ml	85	7.650 ml	+	0.850 ml	8.50 ml
36	3.240 ml	+	0.360 ml	3.60 ml	86	7.740 ml	+	0.860 ml	8.60 ml
37	3.330 ml	+	0.370 ml	3.70 ml	87	7.830 ml	+	0.870 ml	8.70 ml
38	3.420 ml	+	0.380 ml	3.80 ml	88	7.920 ml	+	0.880 ml	8.80 ml
39	3.510 ml	+	0.390 ml	3.90 ml	89	8.010 ml	+	0.890 ml	8.90 ml
40	3.600 ml	+	0.400 ml	4.00 ml	90	8.100 ml	+	0.900 ml	9.00 ml
41	3.690 ml	+	0.410 ml	4.10 ml	91	8.190 ml	+	0.910 ml	9.10 ml
42	3.780 ml	+	0.420 ml	4.20 ml	92	8.280 ml	+	0.920 ml	9.20 ml
43	3.870 ml	+	0.430 ml	4.30 ml	93	8.370 ml	+	0.930 ml	9.30 ml
44	3.960 ml	+	0.440 ml	4.40 ml	94	8.460 ml	+	0.940 ml	9.40 ml
45	4.050 ml	+	0.450 ml	4.50 ml	95	8.550 ml	+	0.950 ml	9.50 ml
46	4.140 ml	+	0.460 ml	4.60 ml	96	8.640 ml	+	0.960 ml	9.60 ml
47	4.230 ml	+	0.470 ml	4.70 ml	97	8.730 ml	+	0.970 ml	9.70 ml
48	4.320 ml	+	0.480 ml	4.80 ml	98	8.820 ml	+	0.980 ml	9.80 ml
49	4.410 ml	+	0.490 ml	4.90 ml	99	8.910 ml	+	0.990 ml	9.90 ml
50	4.500 ml	+	0.500 ml	5.00 ml	100	9.000 ml	+	1.000 ml	10.00 ml

Table 1: 1:10 dilution of conjugate concentrate with Diluent / Wash Buffer Working Dilution

Preparation of the test run using the ViraChip® Software

Basic processes are **Assembling**, **Loading**, **Processing**, **Scanning** and **Analysing**.

After starting the ViraChip® Software you may:

- **Assemble:** Test selection and input of sample data.
- **Load:** Template for preparing the MTP and entering lot specific factors. These are scanned using a 2D bar code scanner from the packaging label of the MTP. **Each MTP can carry only one lot number for each ViraChip® test.**
- **Process:** Data transfer to processor.

Borrelia ViraChip® IgM Test Kit

- 3 -

Assay procedure ^{#)}

1. Place the needed amount of wells into the holding frame. Fill free positions of the last column in the holding frame with empty wells.
2. Add 300 µl Diluent / Wash Buffer Working Dilution to each well and incubate by rocking for approx. 5 minutes, aspirate.
3. Add 100 µl of each diluted patient serum or 100 µl of each diluted control serum.
4. Incubate by rocking for 30 minutes at RT.
5. Aspirate the liquid.
6. 3 x washing:
 - add 300 µl Diluent / Wash Buffer Working Dilution
 - incubate by rocking for 5 minutes at RT
 - aspirate the liquid
7. Add 100 µl Conjugate Working Dilution.
8. Incubate by rocking for 30 minutes at RT.
9. Aspirate the liquid.
10. 3 x washing as in step 6.
11. Add 300 µl distilled or deionised water and incubate by rocking for approx. 5 minutes at RT.
12. Aspirate the liquid.
13. Add 100 µl Chromogen / Substrate Solution.
14. Incubate by rocking for 15 minutes at RT.
15. Stop the reaction by aspirating the liquid.
16. 3 x washing by adding 300 µl distilled or deionised water each time.
17. Dry wells.
18. Measure and interpret wells.

Place the wells into the holding frame accordingly to the layout. Pay attention that no plastic particles fall into the wells while breaking the bars.

Make sure the bottoms of the wells are completely covered with liquid. Use a 3D rocker with a rocking frequency of approx. 750 rpm or a 2D rocker with a rocking frequency of approx. 20 Hz. The aspiration needles must not touch the bottom of the wells.

Add diluted patient sera and diluted control sera directly into the wells.

Make sure the bottoms of the wells are completely covered with liquid. Use a 3D rocker with a rocking frequency of approx. 750 rpm or a 2D rocker with a rocking frequency of approx. 20 Hz. Avoid spilling of liquid.

The aspiration needles must not touch the bottom of the wells.

Make sure the bottoms of the wells are not damaged while adding the Diluent / Wash Buffer Working Dilution. The aspiration needles must not touch the bottom of the wells.

Make sure the bottoms of the wells are completely covered with liquid.

Use a 3D rocker with a rocking frequency of approx. 750 rpm or a 2D rocker with a rocking frequency of approx. 20 Hz.

The aspiration needles must not touch the bottom of the wells.

Make sure the bottoms of the wells are not damaged while adding the Diluent / Wash Buffer Working Dilution. The aspiration needles must not touch the bottom of the wells.

Make sure the bottoms of the wells are completely covered with liquid.

The aspiration needles must not touch the bottom of the wells.

Make sure the bottoms of the wells are completely covered with liquid.

Use a 3D rocker with a rocking frequency of approx. 750 rpm or a 2D rocker with a rocking frequency of approx. 20 Hz.

The aspiration needles must not touch the bottom of the wells.

Without incubation.

Dry the wells under continuous airflow for 20 minutes. Alternatively just by exposure to air for at least 12 hours at RT.

Measurement of spot intensities have to be performed within 24 hours (store MTP in a dark place) e.g. by the CLAIR camera system from Sensovation. The subsequent interpretation is done by the ViraChip® Software.

Assay interpretation with the ViraChip® Software

1. After measuring the spot intensities the interpretation of the ViraChip® Microarrays is performed using the ViraChip® Software. A detailed description of each step can be found in the ViraChip® Software Manual.*
2. Check validity of ViraChip® Microarrays.
The validity check is performed by the ViraChip® Software automatically.
3. Check spot assignment.
The spot layout is shown in Fig. 1. The spot assignment is performed by the ViraChip® Software automatically.

By using the ViraChip® Software you are able to:

- **Scan:** Measurement of the single ViraChip® Microarrays e.g. by the CLAIR camera system (Sensovation).

- **Analyse:** Calculation of the total result from the data

A test run is valid, if the following spots are detectable on each ViraChip® Microarray:

- Serum controls (sc)
- Conjugate controls IgM (ccM)
- Calibrator controls (cal)

and if the following spot is **not** visible:

- Negative control (nc)

If these validation criteria are not fulfilled, the ViraChip® Microarray is classified as invalid. ViraChip® Microarrays that are invalid must not be interpreted and should be repeated.

If multiple conjugate controls are detectable, the strongest spots must indicate the conjugate class being used.

The visual verification of the spots being detected is done by the user. For implausible assignments or wrongly detected spots the QC selection field of the ViraChip® Software has to be changed to "invalid". This sample should be repeated.

^{#)} For automated processing the incubation times of single steps of the procedure may be adjusted to the respective processor type. Refer to section "Notes to Equipment and Software".

^{*)} We gladly provide a ViraChip® Software Manual.

4. Assessment of ViraChip[®] Microarrays.

According to quality laboratory guidelines, the use of a calibrator control or cut off control for each run is recommended (16, 25). The calibrator controls of the Borrelia ViraChip[®] IgM are integrated on the ViraChip[®] Microarray. The assessment is performed by the ViraChip[®] Software automatically.

5. Interpretation of patient spots.

The identified spot triplets of each patient sample have to be considered as symptoms of the disease. A final clinical diagnosis should always be made considering anamnesis, clinical manifestations and laboratory data (24).

The measured mean intensity of the calibrator controls is multiplied by the lot specific factor for each antigen (spot triplet). The resultant value is used as cut-off for the assessment of the respective antigen.

A spot triplet is considered as **distinct** if its mean intensity is equal to or higher than the intensity of the respective cut off. A spot triplet is not assessed if its mean intensity is **lower** than the intensity of the respective cut off or if it is **not present**.

"If the pattern of reactive bands meets the specific conditions, the result is positive, i.e. the positive result of an EIA or another test of the first step is confirmed. If, despite the presence of specific diagnostic bands, the criteria for a positive result are not fulfilled, the result is considered equivocal. In such a case a follow-up control may be recommended." (18)

The following antigens of the Borrelia ViraChip[®] IgM are considered as **highly specific** for Borrelia species for the detection of IgM antibodies: **p41, p39, OspC, p21, Osp17 and VlsE**.

Figure

Antigens:

Each Borrelia specific antigen (**p41, p39, OspC, Osp17 and VlsE**) is spotted three times with the same concentration as spot triplet. Each spot triplet corresponds to one band on an immunoblot.

Controls:

The following integrated controls are implemented on the Borrelia ViraChip[®] IgM:

Serum controls (sc), negative control (nc), conjugate controls (ccG, ccM) and calibrator controls (cal).




Figure1: Schematic drawing of one well of the microtiter plate with the Borrelia ViraChip[®] IgM Microarray (magnified). Spot layout for antigens and integrated controls.

Interpretation criteria

Distinct Borrelia specific spot triplets are calculated in relation to the calibrator control and the lot specific factors by the ViraChip[®] Software. Calibrator Controls are implemented on each ViraChip[®] Microarray.

Identified spot triplet	Result	Interpretation
At least one distinct spot triplets out of: p41, p39, OspC, Osp17 or VlsE	Positive	Specific IgM antibodies against Borrelia species detectable. An infection with Borrelia species is probable.
No or only weak spot triplets	Negative	No specific IgM antibodies against Borrelia species detectable. If an infection is suspected, check a second sample for IgM and IgG specific antibodies after 2-3 weeks.

Rheumatoid factor can affect the reactivity of anigen bands in IgM tests. In the case of unclear band constellation use RF-Absorbent (ViraSorb, 5 ml, Order No.: CB003) for IgM analysis.

Nomenclature and description of Borrelia species bands from literature

Band nomenclature:	Antigen:	Comments:
p41	Limited specificity Flagellum protein	Cross-reactivities with other spirochetes and other flagellum carrying bacteria are described (2,15,22,23).
p39	Highly specific BmpA (Borrelia <u>m</u> embrane protein <u>A</u>)	In many patients antibodies against p39 are already detectable in the early stage of the disease (1,15).
OspC	Highly specific OspC (<u>O</u> uter <u>s</u> urface protein <u>C</u>)	At least 13 different immuno-distinct types of OspC are known. IgM antibodies against OspC are often the first Borrelia specific antibodies detectable in patients and may appear before IgM antibodies against p41 (1,8,12,13,21,22).
Osp17	Specific Osp17 (<u>O</u> uter <u>s</u> urface protein <u>17</u>)	Binding to Decorin on the host cell. Antibodies against Osp17 are described as specific. The presence of IgG antibodies is associated among others with Arthritis and Neuroborreliosis. Species specific (9,10,19).
VlsE	Specific VlsE (Variable major protein (VMP) like sequence Expressed)	Antibodies against VlsE are described as specific. IgM antibodies are already being developed at an early stage and remain until the late stage of the disease (9).

Diagnostic significance of antibodies against Borrelia species

- IgG antibodies** are produced for the first time several weeks to months after infection and are often not detectable in early stages of infection (22). In suspicion of a recent infection, IgM antibodies should be checked and a second sample should be analysed later. Patients in the 2nd or 3rd stage of the disease are usually positive for IgG antibodies. Antibody titres decrease gradually during convalescence (22).
- IgM antibodies** usually appear 2-3 weeks after onset of the disease for the first time (22). Antibody titres often decline several weeks to months after convalescence. But they may also persist up to several years (7,11,20).
- IgA antibodies** are detectable at an early stage of borreliosis in many patients, in some cases earlier than IgM antibodies.
- The immune response and consequently the spot triplet pattern differ from patient to patient. As a general rule: The

number of antibody types and therefore the number of specific spot triplets is increasing with progression of the disease (1).

- An early antibiotic therapy can suppress the development of antibodies (17).

- Medication and immunoglobulin therapy can cause unspecific antibody reactions (24).

- Cross reactivities to Borrelia antigens are described for infections with Treponema, Leptospira and other bacteria with flagella (2,15,22). An acute EBV infection can cause a polyclonal stimulation of Borrelia antibodies (22). If IgM antibodies against OspC or p41 are detected without clinical symptoms for Borreliosis it needs testing for an EBV infection. Cross reactivities in cases of autoimmune diseases, MS, ALS, Influenza and Syphilis are described as well.

Performance data

Sensitivity

106 sera from patients with the Lyme borreliosis at early stages Erythema migrans, Erythema chronicum migrans, multiple Erythema migrans and at late stages Acrodermatitis chronica atrophicans, Lyme arthritis and neuroborreliosis were assessed to determine the sensitivity of the Borrelia ViraChip® IgM Test Kit.

Stage of Lyme borreliosis	Borrelia ViraChip® IgM positive, % (n)	Borrelia ViraChip® IgG / IgM positive/equivocal, % (n)
Early stages (n= 48)	46% (22)	71% (34)
Late stages (n= 58)	47% (55)	97% (56)

To determine the correlation the Borrelia ViraStripe® IgM Test kit was used as a reference test.

Correlation	Borrelia ViraChip® IgM, % (n)
Reference test positive (n= 48)	94% (45)
Reference test negative (n= 58)	93% (54)
Total (n=106)	93% (99)

Borrelia ViraChip® IgM Test Kit

- 6 -

Specificity

To determine the specificity of the Borrelia ViraChip® IgM Test Kit 134 sera from blood donors showing a negative result with a reference test (Borrelia ViraStripe® IgG and IgM Test Kit) were assessed.

Collective	Borrelia ViraChip® IgM negative % (n)	Borrelia ViraChip® IgG / IgM negative % (n)
Blood donors (n= 134)	98% (131)	96% (129)

Warnings and precautions

1. All human serum components were tested for HCV, HIV1,2 antibodies and HBs antigens and found to be negative. Nevertheless, all human kit components as well as the patient samples should be considered as potentially infectious and handled according to safety precautions. While working with potentially infectious/hazardous materials, all national and international rules, regulations, guidelines and laws must be taken into account. This also applies to storage and disposal of chemicals and reagents being used.
2. While working with hazardous or toxic substances/ biological agents precautions have to be applied following national biosafety regulations. Precautions among others are:
 - Do not pipette by mouth.
 - Wear disposable gloves and safety glasses while working.
 - Do not eat, drink or smoke in the working area.

3. The chromogen/substrate solution contains BCIP and NBT. Avoid contact with skin and mucous membranes. In case of contact with skin and eyes wash immediately with large quantities of water.
4. Samples and all potentially contaminated materials must be decontaminated using validated laboratory techniques, e.g. by autoclaving 20 minutes at 121°C. Liquid disposals can be mixed with sodium hypochlorite to a final concentration of 1% sodium hypochlorite. Incubate 30 min for complete disinfection.
5. Please refer to material safety data sheets for detailed information on potential risks, first aid guidelines, accidental release measures, handling and storage recommendations, personal protective equipment, directions for disposal and indications to toxicology.
6. Dust and other contaminations in the wells of the MTP have to be avoided, as this might lead to invalid results.

Storage and stability of reagents

1. **ViraChip® Microarrays:** In closed bags stable until the expiration date if stored at 2-8°C.
2. **Conjugate Concentrate:** Stable until the expiration date if stored at 2-8°C.
3. **Conjugate Working Dilution:** Prepare freshly prior to each run. Do not store for further use.
4. **Diluent / Wash Buffer Concentrate, 10x:** Stable until the expiration date if stored at 2-8°C.

5. **Diluent / Wash Powder:** Stable until the expiration date if stored at 2-8°C.
6. **Diluent / Wash Buffer Working Dilution:** Stable for 2 weeks if stored at 2-8°C. For longer storage, aliquot and freeze at -20°C.
7. **Chromogen / Substrate Solution:** Stable until the expiration date if stored at 2-8°C. Avoid exposure to light!

Specimen indications

1. The **Borrelia ViraChip® IgM Test Kit** must be used with human serum.
2. Specimens must not be microbially contaminated.
3. Using heat-inactivated, haemolytic, icteric or lipemic specimens may lead to invalid results.

4. Normally, human serum can be stored up to 5 days at 2-8°C. Specimens may be stored at -20°C (or below) for long term storage.
5. Prior test processing, specimens should have reached room temperature. Mix specimens carefully after thawing. Precipitates in specimens can be removed by centrifugation.
6. Avoid multiple freeze and thaw cycles.

Limitation of the procedure

1. To ensure reliable results, follow carefully the Instruction for Use and "Good Laboratory Practice".
2. A positive result is based on elevated specific antibody titres and should be considered as a symptom. The correlation to a disease is only conditionally possible.
3. A negative result does not exclude a contact with the pathogen or the presence of a disease.

4. Adequately trained personnel only should perform the assay procedure.
5. The detection of specific antibodies can vary within different assays from different manufacturers and can lead to different results due to different sensitivity, specificity and assay methodologies.
6. ViraChip® Microarrays showing a high background level should not be interpreted, especially if spot intensities are lighter than the background level.

Notes to Equipment and Software

1. Automatized processing requires usage of processor type specific test procedures which are validated and programmed by Viramed Biotech AG.
2. Usage of processor specific consumables requires approval of the respective configurations according to manufacturer's instruction by Viramed Biotech AG.
3. The equipment and software configuration provided by Viramed Biotech AG must not be changed. Any alteration can lead to false results.
4. Only equipment specific software must be used. Changes of configuration data may just be done by Viramed Biotech AG.

5. Only measuring devices approved by Viramed Biotech AG are allowed to be used.
6. Assay interpretation of ViraChip® Microarrays has to be performed using the ViraChip® Software. A manually / visually interpretation is not possible.
7. *In vitro* diagnostics must not be used beyond expiration date as reliable results may not be possible.
8. Efficient washing after each incubation step is essential for consistent results; insufficient washing may lead to false results.

Borrelia ViraChip® IgM Test kit

- 7 -

Literature

1. AGUERO-ROSENFELD, M.: J. Clin. Microbiol. 3090-3095 (1993)
2. ALFEN, I. et al.: Lab. med. 12-19 (1994)
3. CDC/ASTPHLD: Lyme disease Workgroup Recommendations, Dearborn (1994)
4. DITTON H.J.: FEMS Microbiol. 217-230 (1992)
5. DIN 58969-44: Medical microbiology – Diagnostics of infectious diseases in serology and molecular biology, Part 44: Immunoblot (IB); Special requirements for the detection of antibodies against Borrelia burgdorferi (2005)
6. DRESSLER F.: J. Infect. Dis. 392-400 (1993)
7. HAUSER, U. et al.: Interpretation Criteria for Standardized Western Blots for Three European Species of Borrelia burgdorferi Sensu Lato, J. Clin. Microbiol. 1433-1444 (1997)
8. FINGERLE, V. et al.: J. Clin. Microbiol. 1861-1869 (1995)
9. SCHULTE-SPECHTEL, U. et al.: Significant Improvement of the Recombinant Borrelia-specific Immunoglobulin G Immunoblot test by addition of VlsE and a DbpA homologue derived from Borrelia garinii for Diagnosis of Early Neuroborreliosis. J. Clin. Microbiol. (41): 1299-1303 (2003)
10. HEIKKILÄ, T. et al.: Species-Specific Serodiagnosis of Lyme Arthritis and Neuroborreliosis due to Borrelia sensu stricto, B. afzelii, and B. garinii by using Decorin Binding Protein A. J. Clin. Microbiol. (40): 453-460 (2002)
11. HAMMERS-BERGREN, S.: J. Clin. Microbiol. 1519-1525 (1994)
12. JAURIS-HEIPKE, S. et al.: J. Clin. Microbiol. 1860-1866 (1995)
13. JAURIS-HEIPKE, S.: Int. Conference of Lyme Borreliosis. (1994)
14. LAM, T.: Infection and Immunity, 290-298 (1994)
15. MA, B.: J. Clin. Microbiol. 370-376 (1992)
16. RILI-BÄK: Bäk-Richtlinie zur Qualitätssicherung quantitativer laboratoriumsmedizinischer Untersuchungen, (2008), www.bundesärztekammer.de
17. PREAC-MURSIC, V.: Infect. 355-359 (1989)
18. WILSKE, B. et al.: MiQ: Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik 12-2000: Lyme - Borreliose, URBAN&FISCHER, (2000)
19. JAURIS-HEIPKE, S. et al.: Osp 17, a novel immunodominant outer surface protein of Borrelia afzelii: Recombinant expression in Escherichia coli and its use as a diagnostic antigen for serodiagnosis of Lyme borreliosis. Med. Microbiol. Immunol. (Berl), 187 (4): 213-219 (1999)
20. HAUSER, U. et al.: J. Clin. Microbiol. 2241-2247 (1999)
21. WILSKE, B. et al.: Phenotypic Analysis of Outer Surface Protein C (OspC) of Borrelia burgdorferi Sensu Lato by Monoclonal Antibodies: Relationship to Genospecies and OspA Serotype; J. Clin. Microbiol. 103-109 (1995)
22. WILSKE, B.: Diagnose und Labor, (1990)
23. ZÖLLER, L.: J. Clin. Microbiol. 174-182 (1991)
24. THOMAS, L.: Labor und Diagnose, Med. Verlagsgesellschaft Marburg (2008)

Symbols used

	Manufacturer		Order Number
	Refer to Instructions for Use		Use by / Expiration Date
	In-Vitro Diagnostic Medical Device		Temperature Limitation (Storage)
	Lot Number		Positive serum control
	Sufficient for 96 tests		Negative serum control